


# CASE

For Wound Bed Preparation










For better care and wound healing outcomes

# Contents

You can access the relevant section by clicking on one of the links below. Press the  button to return to this page.

	<a href="#">Introduction</a>
	<a href="#">Cause: Identify barriers to wound healing</a>
	<a href="#">Assess: TIMES</a>
	<a href="#">Assess: Wound infection continuum</a>

	<a href="#">Select: T - Tissue</a>
	<a href="#">Select: I - Infection</a>
	<a href="#">Select: M – Moisture</a>
	<a href="#">Select: E – Edge</a>
	<a href="#">Select: S – Surrounding Skin</a>
	<a href="#">Evaluate</a>
	<a href="#">References</a>

## Introduction

The overall goal of wound bed preparation (WBP) is to create an optimal wound healing environment

Wound bed preparation uses the TIMES framework as a structured approach to wound assessment<sup>1</sup>








TIMES identifies barriers to the healing process



Wound bed preparation enhances the effectiveness of therapeutic measures<sup>2,3</sup>

## Cause: Identify barriers to wound healing

Holistic wound assessment considers the 'whole' patient and identifies underlying causes that compromise wound healing








Patient factors to consider:					
	Medical and surgical history		Medication		Activities of daily living
	Nutrition and hydration		Pain assessments		Psychosocial issues
			Overall skin integrity		

## Assess

The five components that underpin Wound Bed Preparation are T.I.M.E.S.

Click on the relevant tab for more information:




	<b>T - Tissue</b>			<b>I - Infection</b>
	<b>S - Surrounding Skin</b>			<b>M - Moisture</b>
	<b>E - Edge</b>			<b>For leg ulcers</b> Adhere to local leg ulcer management policies and guidelines. Carry out a full holistic vascular assessment and consider the need for compression. <sup>1</sup>

**T – Tissue**

Close



- Assess whether there is dead or devitalised tissue in the wound
- This can harbour bacteria and hinder wound healing<sup>4,5</sup>
- Consider whether biofilm may be present, which will need to be managed as part of the debridement process

**T – Tissue****S - Surrounding Skin****I - Infection****M - Moisture****E – Edge****For leg  
ulcers**

Adhere to local leg ulcer management policies and guidelines. Carry out a full holistic vascular assessment and consider the need for compression.<sup>1</sup>

**S - Surrounding Skin**


Close



As well as managing the wound, we need to manage the skin around it. Dry or macerated skin can hinder healing<sup>6</sup>

**Is intervention needed for any of the following?**

- |   |  |  |
|---|--|--|
| <ul style="list-style-type: none"> <li>• Is the skin red or inflamed?</li> </ul>        | <ul style="list-style-type: none"> <li>• Are there any signs of hyperkeratosis or dry skin?</li> </ul> | <ul style="list-style-type: none"> <li>• Are there signs of lipodermasclerosis or excoriation?</li> </ul>    |
| <ul style="list-style-type: none"> <li>• Is there any itching or blistering?</li> </ul> |  | <ul style="list-style-type: none"> <li>• Is there any moisture-associated dermatitis?<sup>4</sup></li> </ul> |

**T - Tissue****S - Surrounding Skin****I - Infection****M - Moisture****E - Edge**For leg  
ulcers

Adhere to local leg ulcer management policies and guidelines. Carry out a full holistic vascular assessment and consider the need for compression.<sup>1</sup>


**E – Edge**

Close



Assessing the edge of the wound can help you to see whether a wound is progressing. The edges should be contracting<sup>6</sup>

Over granulation and rolled edges may be a cause for concern and may require specialist referral

**T – Tissue****S - Surrounding Skin****E – Edge****I - Infection****M - Moisture****For leg  
ulcers**

Adhere to local leg ulcer management policies and guidelines. Carry out a full holistic vascular assessment and consider the need for compression.<sup>1</sup>



## M – Moisture

Close



Exudate is fluid which leaks out of the blood vessels into the surrounding tissue. It is produced throughout the wound healing process

Any exudate needs to be assessed for<sup>1</sup>:



Volume

- High volume can indicate increased bacterial levels
- Low volume can indicate dehydration or ischaemia



Appearance

- e.g. Green exudate can indicate infection (pseudomonas)
- Ideal colour of exudate = white wine



Viscosity


- Thicker viscosity indicates higher protein levels which may indicate infection



T – Tissue



S – Surrounding Skin



I – Infection



M – Moisture



E – Edge

For leg  
ulcers

Adhere to local leg ulcer management policies and guidelines. Carry out a full holistic vascular assessment and consider the need for compression.<sup>1</sup>

## I - Infection

Close



Assess the wound for signs and symptoms of infection or inflammation

## Consider the below criteria for the presence of biofilm

[Click here for Biofilm Definition](#)

Failure of appropriate antibiotic treatment

Low-level chronic inflammation

Response to antimicrobial treatment ineffective

Low-level erythema

Recurrence of delayed healing on cessation of antibiotic treatment

Poor granulation/friable hypergranulation

Delayed healing despite optimal wound management and health support

Secondary signs of infection


Increased exudate/moisture

[Click here for Infection Continuum](#)


T - Tissue



S - Surrounding Skin



I - Infection



M - Moisture



E - Edge

For leg ulcers


Adhere to local leg ulcer management policies and guidelines. Carry out a full holistic vascular assessment and consider the need for compression.<sup>1</sup>

## Biofilm definition

Close



Biofilms are bacteria which attach to surfaces. Biofilms produce their own encapsulation that is tolerant to antimicrobial agents (this includes antibiotic and antimicrobial)




World Union of Wound Healing Societies (WUWHS), Florence Congress, Position Document. Management of Biofilm. Wounds International 2016

## Assess: I - Infection

The **wound infection continuum** provides a framework which looks at the impact microbes have on a wound and wound healing.

Please select the options in pink for more info



International Wound Infection Institute (IWII) (2016) Wound infection in clinical practice. Wounds International, London. Available at: <https://www.woundsinternational.com/resources/details/iwiiwound-infection-clinical-practice> (accessed 24.03.20)

## Assess: I - Infection



### Contamination

Close



- All wounds may acquire micro-organisms. Host defenses will stop them multiplying
- At this stage - their presence is only transient and wound healing is not delayed

Contamination

Colonisation

Local infection

Spreading infection

Systemic infection

←--- Vigilance required ---→ ←--- Intervention required ---→

No antimicrobials indicated


Topical antimicrobial

Systemic and topical antimicrobials

[< Back to Infection Page](#)



International Wound Infection Institute (IWII) (2016) Wound infection in clinical practice. Wounds International, London. Available at: <https://www.woundsinternational.com/resources/details/iwiiwound-infection-clinical-practice> (accessed 24/03/20)

## Assess: I - Infection



International Wound Infection Institute (IWII) (2016) Wound infection in clinical practice. Wounds International, London. Available at: <https://www.woundsinternational.com/resources/details/iwiiwound-infection-clinical-practice> (accessed 24/03/20)

## Assess: I - Infection


**Local infection**
Close 

Covert (subtle signs of local infection)		
Hypergranulation (excessive 'vascular' tissue)	Epithelial bridging and pocketing in granulation tissue	Delayed wound healing beyond expectations
Bleeding, friable granulation tissue	Wound breakdown and enlargement	New or increasing pain
Increasing malodour		

Overt (classic signs of local infection)		
Erythema	Local warmth	Swelling
Purulent discharge	Delayed wound healing beyond expectations	New or increasing pain
Increasing malodour		

Contamination

Colonisation

Local infection

Spreading infection

Systemic infection

←
Vigilance required
→
←
Intervention required
→

No antimicrobials indicated

Topical antimicrobial

Systemic and topical antimicrobials

[< Back to Infection Page](#)

International Wound Infection Institute (IWII) (2016) Wound infection in clinical practice. Wounds International, London. Available at: <https://www.woundsinternational.com/resources/details/iwiiwound-infection-clinical-practice> (accessed 24/03/20)

## Assess: I - Infection



### Spreading infection

Close



- Extending in duration + / - erythema
- Lymphangitis
- Crepitus
- Wound breakdown / dehiscence with or without satellite lesions
- Malaise / lethargy or non-specific general deterioration
- Loss of appetite
- Inflammation / swelling of lymph glands

Contamination

Colonisation

Local infection

Spreading infection

Systemic infection

←--- Vigilance required ---→ ←--- Intervention required ---→

No antimicrobials indicated

Topical antimicrobial


Systemic and topical antimicrobials

[< Back to Infection Page](#)

International Wound Infection Institute (IWII) (2016) Wound infection in clinical practice. Wounds International, London. Available at: <https://www.woundsinternational.com/resources/details/iwiiwound-infection-clinical-practice> (accessed 24/03/20)



## Assess: I - Infection





International Wound Infection Institute (IWII) (2016) Wound infection in clinical practice. Wounds International, London. Available at: <https://www.woundsinternational.com/resources/details/iwiiwound-infection-clinical-practice> (accessed 24.03.20)

## Select: T - Tissue

If dead or devitalised tissue has been identified through the wound assessment then you need to consider debridement within your management plan

### Aims of debridement<sup>7</sup>:

 Remove:
<ul style="list-style-type: none"> <li>• Necrotic, devitalised, sloughy tissue</li> </ul>
<ul style="list-style-type: none"> <li>• Sources of infection, inflammation</li> </ul>
<ul style="list-style-type: none"> <li>• Exudate, dried exudate and dry skin/hyperkeratosis</li> </ul>
<ul style="list-style-type: none"> <li>• Pus</li> </ul>
<ul style="list-style-type: none"> <li>• Haematoma</li> </ul>
<ul style="list-style-type: none"> <li>• Debris or foreign bodies</li> </ul>
<ul style="list-style-type: none"> <li>• Any other barriers to healing</li> </ul>

 Decrease:
<ul style="list-style-type: none"> <li>• Odour</li> </ul>
<ul style="list-style-type: none"> <li>• Excess moisture</li> </ul>
<ul style="list-style-type: none"> <li>• Risk of infection</li> </ul>

 Stimulate:
<ul style="list-style-type: none"> <li>• Wound edges and epithelialisation</li> </ul>

 Improve:
<ul style="list-style-type: none"> <li>• Quality of life</li> </ul>

**NB.** It is important to know when to refer to a specialist best qualified to debride (e.g DFU). Remember that NOT debriding / referring can potentially cause harm to your patients; involve your multidisciplinary team where you need to in order to provide the care your patients need<sup>6</sup>.

## Select: T - Tissue



### Debridement

Debridement methods require varying levels of expertise. You need to consider your skills and competency to perform the task and refer to a specialist if necessary<sup>5</sup>.









### Debridement Methods

Select the appropriate debridement methods










	Autolytic
	Mechanical
	Larval
	Sharp
	Surgical
	Ultrasonic
	Hydrosurgical








## Select: T - Tissue

<b>Autolytic</b>	<b>Close</b> 		<b>Autolytic</b>
<ul style="list-style-type: none"><li>• Using moisture to soften hard necrotic tissue and liquify slough e.g. hydrogels, hydrocolloids, hydrofibres, but the process may be slow</li><li>• Suitable for self-care by patients and carers</li></ul>			<b>Mechanical</b>
			<b>Larval</b>
			<b>Sharp</b>
			<b>Surgical</b>
			<b>Ultrasonic</b>
			<b>Hydrosurgical</b>








## Select: T - Tissue

<b>Mechanical</b>	Close X		Autolytic
<ul style="list-style-type: none"><li>• Use monofilament pads or debridement cloths containing a surfactant. Monofilament pads can help remove bacteria and biofilm</li><li>• Suitable for self-care by patients and carers</li></ul>		Mechanical	
		Larval	
		Sharp	
		Surgical	
		Ultrasonic	
		Hydrosurgical	








## Select: T - Tissue

<b>Larval</b>	Close X		Autolytic
<ul style="list-style-type: none"><li>• Sometimes referred to as biosurgical and uses the sterile larvae of the greenbottle fly</li><li>• Rapid debridement, but not always acceptable with patients or healthcare workers</li></ul>			Mechanical
			<b>Larval</b>
			Sharp
			Surgical
			Ultrasonic
			Hydrosurgical

## Select: T - Tissue








<b>Sharp</b>	<b>Close</b> ✕		<b>Autolytic</b>
<ul style="list-style-type: none"><li>• Using a scalpel, scissors or curette to remove tissue. Often used in conjunction with other methods</li><li>• Requires specialist assessment prior to implementation</li></ul>			<b>Mechanical</b>
			<b>Larval</b>
			<b>Sharp</b>
			<b>Surgical</b>
			<b>Ultrasonic</b>
			<b>Hydrosurgical</b>

## Select: T - Tissue








<b>Surgical</b>	Close X		Autolytic
<ul style="list-style-type: none"><li>Excision and removal performed by a specialist in an acute, operating theatre environment</li></ul>		Mechanical	
		Larval	
		Sharp	
		Surgical	
		Ultrasonic	
		Hydrosurgical	



## Select: T - Tissue

<b>Ultrasonic</b>	Close X
<ul style="list-style-type: none"><li>• Using ultrasound either directly on the wound bed or by an atomised solution</li><li>• Requires specialist assessment and implementation</li></ul>	
	 Autolytic
	 Mechanical
	 Larval
	 Sharp
	 Surgical
	 Ultrasonic
	 Hydrosurgical

## Select: T - Tissue

<b>Hydrosurgical</b>	Close X		Autolytic
<ul style="list-style-type: none"><li>• The use of a high energy saline beam as a cutting implement</li><li>• Requires specialist assessment and implementation</li></ul>		Mechanical	
		Larval	
		Sharp	
		Surgical	
		Ultrasonic	
		Hydrosurgical	

## Select: I - Infection

### Biofilm management

Removing biofilm is part of the wound bed preparation process. This is an essential step to facilitate healing<sup>5</sup>

A pathway for the management of biofilm:



Strong cleaning efficiency<sup>8,9</sup>



Effective bacterial binding<sup>10</sup>



>99% biofilm removal<sup>9</sup>

Debridement, using monofilament pads, can help remove bacteria and biofilm.

Static chronic wound



Suspected Biofilm



Reduce biofilm - debridement/vigorous cleansing



Prevent recontamination - topical antimicrobials  
suppress biofilm reformation - repeat debridement




Reassess

(Phillips et al, 2010)

## Select: I - Infection

The **wound infection continuum** provides a framework which looks at the impact microbes have on a wound and wound healing.

Managing wound infection: **Please select the options in pink for more info**




International Wound Infection Institute (IWII) (2016) Wound infection in clinical practice. Wounds International, London. Available at: <https://www.woundsinternational.com/resources/details/iwiiwound-infection-clinical-practice> (accessed 24/03/20)

## Select: I - Infection

The **wound infection continuum** provides a framework which looks at the impact microbes have on a wound and wound healing.

Managing wound infection: **Please select the options in pink for more info**




International Wound Infection Institute (IWII) (2016) Wound infection in clinical practice. Wounds International, London. Available at: <https://www.woundsinternational.com/resources/details/iwiiwound-infection-clinical-practice> (accessed 24/03/20)

## Select: I - Infection

The **wound infection continuum** provides a framework which looks at the impact microbes have on a wound and wound healing.





Managing wound infection: **Please select the options in pink for more info**



International Wound Infection Institute (IWII) (2016) Wound infection in clinical practice. Wounds International, London. Available at: <https://www.woundsinternational.com/resources/details/iwiiwound-infection-clinical-practice> (accessed 24/03/20)

## Select: M – Moisture:

When selecting the appropriate product the following criteria should be considered

	<b>Fluid handling capacity</b>	Select the appropriate dressing suitable for the level of exudate
	<b>Correct size and shape</b>	Select the dressing appropriate for the anatomical location and size of the wound
	<b>Dressing fixation</b>	<p>Consider whether you will use bordered dressings or a secondary fixation such as bandages.</p> <p>The dressing should be secure yet provide atraumatic dressing removal</p>
	<b>Retention of fluid under compression</b>	Select the appropriate dressing to handle the level of exudate under compression





## Select: M – Moisture:

Effective exudate management can promote healing, improve quality of life and enhance healthcare effectiveness. Absorbent products vary in the materials they are made from and in their ability to manage exudate. Knowing how they manage fluid is key to selecting the most appropriate and effective dressing/technology for each wound.

### Exudate management dressings

Click to find out about the different exudate management dressings available



	<b>Foams</b>
	<b>Superabsorbers</b>
	<b>Gelling fibres</b>
	<b>Hydrocolloids</b>







## Select: M – Moisture:

Foams		Close X
<b>Composition:</b> <ul style="list-style-type: none"><li>• Vary in thickness</li><li>• With or without a silicone wound contact layer</li><li>• Bordered or non-bordered options</li></ul> <b>Action:</b> <ul style="list-style-type: none"><li>• Absorb exudate, allowing evaporation to occur via a polyurethane top film</li></ul>	<b>Advantages:</b> <ul style="list-style-type: none"><li>• Easy to apply</li><li>• Minimises trauma and pain (for silicone wound contact layer products)</li><li>• Suitable for many wound types</li></ul> <b>Wound characteristics:</b> <ul style="list-style-type: none"><li>• Traditionally, thinner foams have been designed for lower exudate levels</li><li>• More absorbent (thicker) foams can be used for highly exuding wounds</li></ul>	

	Foams
	Superabsorbers
	Gelling fibres
	Hydrocolloids

## Select: M – Moisture:

Superabsorbers		Close X
<b>Composition:</b> <ul style="list-style-type: none"><li>• Multi-layered polymer construction</li></ul> <b>Action:</b> <ul style="list-style-type: none"><li>• Wick moisture from the wound and lock fluid inside the dressing</li></ul>	<b>Advantages:</b> <ul style="list-style-type: none"><li>• Enhanced absorbency</li><li>• Longer wear times</li><li>• Less-frequent dressing changes</li></ul> <b>Wound characteristics:</b> <ul style="list-style-type: none"><li>• Heavily exuding wounds</li></ul>	

	Foams
	Superabsorbers
	Gelling fibres
	Hydrocolloids





## Select: M – Moisture:

Gelling fibres/Alginates		Close	X
<b>Composition:</b> <ul style="list-style-type: none"><li>• 100% carboxymethylcellulose (CMC)</li><li>• 100% alginate</li><li>• Or a combination</li></ul> <b>Action:</b> <ul style="list-style-type: none"><li>• Transforms into a moist, gel-like sheet or conformable gel when absorbing exudate</li><li>• Transmits water from the wound surface</li></ul>	<b>Advantages:</b> <ul style="list-style-type: none"><li>• Maintain moist wound environment</li><li>• Comfortable</li><li>• Conforms to wound</li><li>• Can be used in deep wounds</li></ul> <b>Wound characteristics:</b> <ul style="list-style-type: none"><li>• Moderate to heavily exuding wounds</li><li>• Not to be used on fragile skin</li></ul>		

	<b>Foams</b>
	<b>Superabsorbers</b>
	<b>Gelling fibres</b>
	<b>Hydrocolloids</b>

## Select: M – Moisture:

Hydrocolloids		Close X
<b>Composition:</b> <ul style="list-style-type: none"><li>• Gel-forming agents</li><li>• Waterproof backing, usually polyurethane</li><li>• Available in various shapes and thicknesses</li><li>• Bordered or non-bordered</li></ul> <b>Action:</b> <ul style="list-style-type: none"><li>• Absorbs exudate and forms a gel</li></ul>	<b>Advantages:</b> <ul style="list-style-type: none"><li>• Maintains a moist wound environment</li><li>• Stimulates granulation</li><li>• Does not stick to the wound</li><li>• Are flexible to mould around skin and body contours</li></ul> <b>Wound characteristics:</b> <ul style="list-style-type: none"><li>• Light to moderately exuding wounds</li></ul>	

	<b>Foams</b>
	<b>Superabsorbers</b>
	<b>Gelling fibres</b>
	<b>Hydrocolloids</b>

## Select: E – Edge

Where wound edges are rolled management of the wound is going to be complex and may need early specialist referral

Indicators of non-advancement:

Rolled edges




Crusting

## Select: S – Surrounding Skin

Surrounding skin may be compromised by the current wound management plan

There may be signs of:

☞ Maceration, excoriation, oedema



☞ Dry skin, hyperkeratosis



### For leg ulcers

Adhere to local leg ulcer management policies and guidelines. Carry out a full holistic vascular assessment and consider the need for compression

## Select: S – Surrounding Skin

Surrounding skin may be compromised by the current wound management plan

There may be signs of:

👉 Maceration, excoriation, oedema

Close X

**Maceration, excoriation and oedema  
can be prevented by effective exudate  
management.**

👉 Dry skin, hyperkeratosis



### For leg ulcers


Adhere to local leg ulcer management policies and guidelines. Carry out a full holistic vascular assessment and consider the need for compression

## Select: S – Surrounding Skin

Surrounding skin may be compromised by the current wound management plan

There may be signs of:

👉 Maceration, excoriation, oedema



👉 Dry skin, hyperkeratosis

Close X

**Manage and prevent hyperkeratosis through good skin care.**

**Washing, drying and application of emollients should be an integrated part of every patients care plan**

**For leg ulcers**

Adhere to local leg ulcer management policies and guidelines. Carry out a full holistic vascular assessment and consider the need for compression



## Evaluate

Wound healing is a dynamic process, meaning that your assessment and management needs to be constantly evaluated. If current treatment plan is not working, consider a reassessment of TIMES and modify the care plan as needed.

### Supported self care

On consideration of a patients ability to participate in supported self care refer to CASE - Assessing a patients ability to self care.

Patients need to know what to look for and when to ask for help

These might include:

- Unusual wound leakage, pain or smell
- Wound increasing in size
- Increasing redness around the wound
- Generally feeling unwell

Patients can be advised to photograph their wound to measure its progress.



## References

1. Wounds UK (2016) Best Practice Statement: Holistic Management of Venous Leg Ulceration. Wounds UK, London  
Available online: [www.wounds-uk.com](http://www.wounds-uk.com)
2. Atkins L, Bucko Z, Conde Montero E et al (2019) Implementing TIMERS: The race against hard-to-heal wounds. *Journal of Wound Care* 28(Suppl 3) S1-S49
3. Falanga V. Classifications for wound bed preparation and stimulation of chronic wounds. *Wound Repair Regen.* 2000; 8(5): 347-52.
4. Atkin L (2016) A guide to wound debridement. Independent Nurse. MA healthcare Ltd, London. 04 February, 1-4
5. Vowden K and Vowden P (2011) Debridement made easy. *Wounds UK* 7(4) S1-S4
6. Harries R, Bosanquet D, Harding K (2016) Wound bed preparation: TIME for an update. *International Wound Journal* 13(Suppl.S3) 8-14
7. Strohal R, Apelqvist J, Dissemond J, et al (2013) EWMA document: Debridement. *J Wound Care* 22(Suppl 1): S1-S52
8. JB4 40031908, In vitro evaluation of removal of viable biofilm by Cutimed® Debriclean and Cutisoft Cotton. August 2018, BSN medical data on file.
9. JB4 40031810, In vitro evaluation of the cleansing effect of Cutimed® Debriclean and Cutisoft Cotton. August 2018, BSN medical data on file.
10. JB4 40031811, In vitro evaluation of the bacterial adhesion to Cutimed® Debriclean and Cutisoft Cotton. July 2018, BSN medical data on file.
11. NHS Clinical Evaluation Team (2018): Clinical Review: Superabsorbent Dressings.