

JOBST® Comprifore

Multi-layer compression bandage system

To suit the needs of both clinician and patient, a compression bandage system has to provide effective treatment whilst being easy to use. JOBST Comprifore is applied in a simple technique to deliver safe and sustained compression.

Application



After applying a suitable wound dressing, start with JOBST Comprifore #1 at the base of the toes and anchor the bandage with 2 turns. Apply the padding as smoothly and evenly as possible without stretching.



Ensure the heel is covered and work up the leg in a spiral technique with a 50% overlap. Finish below the knee.



Proceed in just the same way with JOBST Comprifore #2: Ensure the bandage is firm at the arch of the foot.



If using JOBST Comprifore lite for reduced compression, go to picture 5. Apply JOBST Comprifore #3 in figure of 8 technique with 50% stretch using the coloured central line as guidance to achieve the required 50% overlap.



After giving anchor to JOBST Comprifore #4, go round the ankle and work up the leg in a spiral technique with a 50% overlap at 50% stretch.



The cohesive material will form a permanent bond and secure the bandage. Smooth out the layers with your hands and examine for any gaps.