



Animal Healthcare Catalogue



About Essity

We are focused on the development of world-class, branded products that offer high quality solutions.

Essity is a global leader in the worldwide healthcare market, specialising in the areas of Wound Care and Fracture Management.

Essity stands for brands such as **Soffban®**, **Tensoplast®**, **Co-Plus®**, **Gypsona®** and **Easifix®** with a long history that you have come to trust.

Whether it is dealing with an acute fracture or simply ensuring that an animal is bandaged effectively, Essity has the products to suit all of your veterinary healthcare needs.

Please visit our website for further product information - www.medical.essity.co.uk

Key Account Contacts

Please contact **Charlotte Wood** and **Francesca Gordon** for more information

Tel: **01482 670124**

Email: animalhealthcare@essity.com



Introducing innovative advanced wound care to veterinary professionals



Cutimed® Siltec®



Cutimed® Sorbion®



Leukomed® Sorbact®



Cutimed® Gel

Contents

About Essity	2	Bandages	
Contact Information	2	Padding	
		Soffban Synthetic	12
		Soffban Natural	12
		Retention	
		Easifix	12
		Easifix K	12
		Elastic Adhesive	
		Tensoplast	13
		Cohesive	
		Co-Plus	14
		Easifix Cohesive	14
		Support	
		Delta Terry-Net C	15
Wound dressings		Tapes	
Exudate management		Hypafix	18
Cutimed Siltec	5	Hypafix Transparent	18
Cutimed Sorbion	6	Leukopor	18
Infection management		Leukosilk	18
Cutimed Sorbact	7	Leukofix	19
Wound bed management		Leukoplast Sleek	19
Cutimed Debriclean	8	Tensoplast	19
Cutimed Gel	8	Strappal	19
Tissue Management			
Cutimed Epiona	9		
Other			
Leukomed Sorbact	10		
Leukomed Control	10		
Leukomed	11		
Leukomed T	11		
Leukomed T plus	11		
Cuticell Classic	11		
Cuticell Contact	11		
Casting Materials		Swabs	
Delta-Cast Elite	16	Propax	20
Delta-Cast Conformable	16	Multisorb	20
Gypsona	17	Fivex	20
Dynacast Prelude	17		
		A Quick Summary Guide	22

What is Cutimed®?

Cutimed is the Essity brand of products designed to facilitate the wound healing process



Cutimed® DebricClean

Mechanical debridement pad



Cutimed® Siltec®

Foam dressing for low to moderate exudate levels



Cutimed® Gel

Hydrogel for debriding necrotic or sloughy tissue



Cutimed® Sorbion®

Superabsorbent pad for moderate to high exudate levels



Cutimed® Sorbact®

Safe and effective removal of bacteria to reduce the bioburden in wounds and facilitate the healing process



Cutimed® Epiona

Collagen dressing to help boost chronic wounds

Cutimed® Siltec®

A range of foam dressings designed to manage any low to moderately exuding wounds

The combination of a silicone wound contact layer and superabsorbers above the foam core mean that Cutimed Siltec is excellent in its exudate management and minimises wound bed damage when removing.

The range

Secure, comfortable, conformable

Cutimed Siltec products are designed for atraumatic dressing changes with minimised pain and to protect against odour. Additionally, the border products are fully water resistant.



Cutimed® Siltec® B

- Fast application, secondary fixation is not necessary



Cutimed® Siltec® Plus

- Strong initial tack for that helping hand effect, easy to apply and easy to re-adjust



Cutimed® Siltec®

- Skin protection with a feather tack
- Excellent exudate management



Cutimed® Siltec® L

- Soft and conformable wound dressing
- Protective dressing for compromised and fragile skin



	Code	Dressing size (cm)	Pack Qty
Cutimed Siltec B (Border) ☒☒☒☒☒☒	73284-00	7.5 x 7.5	1 x 10
	73284-01	12.5 x 12.5	1 x 10
	73284-02	15 x 15	1 x 10
	73284-03	17.5 x 17.5	1 x 5
	73284-04	22.5 x 22.5	1 x 5
	73284-05	10 x 10	1 x 10
	73284-06	10 x 22.5	1 x 10
Cutimed Siltec B Oval ☒☒☒☒☒☒	73284-14	7 x 10 / oval	1 x 10
Cutimed Siltec Plus ☒☒☒☒☒☒	73288-00	5 x 6	1 x 10
	73288-01	10 x 10	1 x 10
	73288-02	10 x 20	1 x 10
	73288-03	15 x 15	1 x 10
	73288-04	20 x 20	1 x 5
Cutimed Siltec ☒☒☒☒☒☒	73285-00	5 x 6	1 x 10
	73285-01	10 x 10	1 x 10
	73285-02	10 x 20	1 x 10
	73285-03	15 x 15	1 x 10
	73285-04	20 x 20	1 x 5
Cutimed Siltec L (Light) ☒☒☒☒☒☒	73283-00	5 x 6	1 x 10
	73283-01	10 x 10	1 x 10
	73283-02	15 x 15	1 x 10

Wound depth

- ☒ Superficial
- ☒ Deep
- ☒ Superficial + deep

Wound phase

- ☒ Necrotic
- ☒ Infected
- ☒ Sloughy
- ☒ Granulating
- ☒ Epithelializing

Exudate level

- ☒ Dry to low
- ☒ Low to moderate
- ☒ Moderate to high
- ☒ High to excessive

Cutimed® Sorbion®

Super-absorbent dressings

Highly to excessively exuding wounds are hard to be controlled - even superior foam dressings such as Cutimed Siltec need to be changed quite often for those wounds

Consequences of inappropriate dressings are;

- Excessive wound exudate can cause skin maceration around a wound, which may delay healing and lead to other complications.¹ A cycle of moisture, friction, bacteria and breakdown may then be established²
- Frequent dressing changes due to low capacity could damage the wound bed³

Cutimed® Sorbion® has demonstrated its efficacy in the treatment of highly exuding wounds. Utilising patented Hydration Response Technology, Cutimed® Sorbion® has set standards in modern exudate management.



Expansion rim

- Broad wound contact thanks to even expansion, even under high saturation

Ultrasonic sealing

- Elimination of adhesives and binders for a particularly well-tolerated wound care - mean less patient allergies

Polymer gelling

- Moist wound environment
- Ensures close contact with the wound bed

Fibre matrix


- Avoids clumping
- Vertical absorption, safe and fast

Fibre matrix with polymers

- High exudate absorption capacity
- Extension of dressing change intervals
- No rewet into the wound
- Reduces the risk of maceration

Polypropylene cover

- Strong tear resistance
- High breathability



Code	Size	Pack Qty	Case Qty
Cutimed Sorbion Sachet EXTRA			
73234-00	5cm x 5cm	1 x 5	115 x 5
73234-01	7.5cm x 7.5cm	1 x 5	76 x 5
73234-02	10cm x 10cm	1 x 5	44 x 5
73234-03	20cm x 10cm	1 x 5	35 x 5
73234-04	20cm x 20cm	1 x 5	21 x 5
73234-05	30cm x 20cm	1 x 5	12 x 5
Cutimed Sorbion Sachet XL			
73240-01	45cm x 25cm	1 x 5	16 x 5



Wound depth

- Superficial
- Deep
- Superficial + deep

Wound phase

- Necrotic
- Infected
- Sloughy
- Granulating
- Epithelializing

Exudate level

- Dry to low
- Low to moderate
- Moderate to high
- High to excessive

1) Cutting and White, 2002;
 2) Fiers, 1996;
 3) European Pressure Ulcer Advisory Panel's Pressure Ulcer Treatment Guidelines;

Cutimed® Sorbact®

Bacteria binding dressings coated with Sorbact® technology

Safely remove bacteria -

Killing the bacteria releases endotoxins from dead cells and cell debris is left in the wound. So unlike antimicrobial dressings that actively kill bacteria, Cutimed® Sorbact® removes it, irreversibly binding it to its surface to reduce bioburden and support wound healing. Cutimed® Sorbact® has no contraindications and a low risk of allergies.

Cutimed® Sorbact® is versatile and effective in a variety of wound stages. It can be used to safely reduce bioburden in the short term, and also in the long-term treatment of already infected wounds.

Cutimed Sorbact dressing pad

An absorbent dressing pad covered with a Sorbact technology coated acetate swab. For wounds with low to moderate exudate levels



Cutimed Sorbact Swabs

Sorbact technology coated acetate wound contact layer. Primary wound dressing for either deep or superficial wounds



Cutimed Sorbact Ribbon

Made from Sorbact technology coated cotton. For use in deep or cavity wounds or wrapping around digits



Cutimed Sorbact Gel

A Cutimed Sorbact swab coated with an amorphous hydrogel for dry, sloughy or low exuding wounds



Cutimed Sorbact Hydroactive

A unique combination of a Cutimed Sorbact swab combined with an absorbent hydropolymer gel matrix. These dressings reduce bacterial load, absorb and lock in wound exudate whilst maintaining a moist wound environment. For low to moderate exudate levels



Cutimed Siltec Sorbact

Cutimed Siltec Sorbact combines infection management with the excellent fluid handling of Cutimed Siltec. For wounds with moderate to high, viscous exudate



Code	Size	Pack Qty	Case Qty
Cutimed Sorbact dressing pads			
72161-01	7cm x 9cm	1 x 5	16 x 5
72162-01	10cm x 10cm	1 x 5	16 x 5
72163-00	10cm x 20cm	1 x 20	2 x 20
Cutimed Sorbact Swabs*			
72164-01	4cm x 6cm folded (11cm x 16cm unfolded)	1 x 5	24 x 5
72165-01	7cm x 9cm folded (17cm x 27cm unfolded)	1 x 5	24 x 5
Cutimed Sorbact Ribbon			
72166-00	2cm x 50cm	1 x 20	6 x 20
72167-00	5cm x 200cm	1 x 10	6 x 10
Cutimed Sorbact Hydroactive (wound pad)			
72646-03	7cm x 8.5cm (5cm x 6.5cm)	1 x 10	10 x 10
72646-04	14cm x 14cm (10cm x 10cm)	1 x 10	10 x 10
72646-05	19cm x 19cm (15cm x 15cm)	1 x 10	10 x 10
72646-06	24cm x 24cm (20cm x 20cm)	1 x 5	10 x 5
72646-07	14cm x 24cm (10cm x 20cm)	1 x 10	10 x 10
Cutimed Sorbact Hydroactive B (wound pad)			
79933-00	7cm x 8.5cm (5cm x 6.5cm)	1 x 10	10 x 10
79933-02	14cm x 14cm (10cm x 10cm)	1 x 10	10 x 10
79933-03	19cm x 19cm (15cm x 15cm)	1 x 10	10 x 10
79933-04	14cm x 24cm (10cm x 20cm)	1 x 10	10 x 10
Cutimed Siltec Sorbact (wound pad)			
79929-00	7.5cm x 7.5cm (3.9cm x 3.9cm)	1 x 10	10 x 10
79929-01	12.5cm x 12.5cm (8.2cm x 8.2cm)	1 x 10	10 x 10
79929-02	15cm x 15cm (10cm x 10cm)	1 x 10	10 x 10
79929-03	17.5cm x 17.5cm (12.3cm x 12.3cm)	1 x 5	10 x 5
79929-04	22.5cm x 22.5cm (16.5cm x 16.5cm)	1 x 5	10 x 5
79929-05	Sacrum 17.5cm x 17.5cm	1 x 5	10 x 5
79929-06	Sacrum 23cm x 23cm	1 x 5	10 x 5
Cutimed Sorbact Gel			
72611-00	7.5cm x 7.5cm	1 x 10	36 x 10
72611-01	7.5cm x 15cm	1 x 10	20 x 10

Wound depth

- Superficial
- Deep
- Superficial + deep

Wound phase

- Necrotic
- Infected
- Sloughy
- Granulating
- Epithelializing

Exudate level

- Dry to low
- Low to moderate
- Moderate to high
- High to excessive

Cutimed® DebriClean

Fast and effective wound debridement

In chronic wounds, devitalised tissue and biofilm forms a physical barrier that inhibits the formation of new tissue. For the wound to heal properly, the devitalised tissue and biofilm must be removed.¹ That's why every layer of the Cutimed® DebriClean pad has been specially designed to support you in the mechanical debridement process and lead to more effective wound healing.



Cutimed® DebriClean

Mechanical Debridement Pad

- A dual action debridement pad with soft and abrasive fibres
- Effectively cleans wounds of slough, biofilms and microorganisms
- Easy and comfortable to use

Code	Size	Pack Qty
76170-00	9cm x 13cm	1 x 5



Cutimed® Gel

Amorphous hydrogel

- A clear, amorphous hydrogel which can be used to debride necrotic and sloughy tissue
- Easy to apply with the sterile nozzle provided, a spatula or directly from the tube
- Compatible with any type of secondary dressing, dependent upon exudate levels

Ideal for: autolytic debridement, sloughy or necrotic tissue, and moisture donation for dry wounds

Code	Size	Pack Qty	Case Qty
72610-00	8g	1 x 10	8 x 10
72610-01	15g	1 x 10	8 x 10
72610-02	25g	1 x 10	8 x 10




Cutimed® Epiona

Boosts healing in chronic and stagnating wounds

Natural mode of action - made from 90% native collagen with 10% alginate

Resorbable - no need to remove from the wound

Reduces wound size by up to 40% in 3 weeks¹



- Accelerated wound healing¹
- Safe application^{1,2}
- Easy to combine with your current therapies (e.g. compression therapy)¹
- Excellent hemostatic properties^{3,4}
- Broad indication for treatment of all wounds healing by secondary intention, free of necrotic tissue^{1,2}

Code	Size	Pack Qty
73229-00	5 x 5cm	1 x 10
73229-01	10 x 10cm	1 x 10



1. Sabo, M, et al., 'A Post-marketing Surveillance Study of Chronic Wounds Treated With a Native Collagen Calcium Alginate Dressing', Ostomy Wound Management, 2018
 2. Cutimed® Epiona - Case Report, Data on File.
 3. Wiegand C., et al., 'A novel native collagen dressing with advantageous properties to promote physiological wound healing', Journal of Wound Care, 2016.
 4. Jesty, J. et al., 'Assessment In Vitro of the Active Hemostatic Properties of Wound Dressings', Journal of Biomedical Materials Research Part B: Applied Biomaterials, 2008

Leukomed® Sorbact®

Innovative post-operative dressing to help to reduce wound bacteria colonisation by its physical mode of action. For post-operative and traumatic wounds

Unique bacteria binding Sorbact®

- Effectively binds hydrophobic bacteria and fungi
- Development of bacterial and fungal resistance is not expected

Bacteria proof adhesive film

- Effectively protects against external contamination
- Breathable and shower-proof



Code	Size	Pack Qty	Case Qty
76199-00	5cm x 7.2cm	1 x 20	56 x 20
76199-01	8cm x 10cm	1 x 20	12 x 20
76199-03	8cm x 15cm	1 x 20	8 x 20
76199-04	10cm x 20cm	1 x 20	12 x 20
76199-05	10cm x 25cm	1 x 20	10 x 20
76199-06	10cm x 30cm	1 x 20	8 x 20
76199-07	10cm x 35cm	1 x 20	8 x 20

Leukomed® Control

Unique transparent post-operative dressing with absorbent hydropolymer pad

- Ultra thin and conformable film for a second skin feeling
- Absorbent hydropolymer gel wound pad
- Completely transparent composition
 - Allows for early detection of wound infection

Ideal for: surgical incisions, lacerations, cuts and abrasions



Code	Size	Pack Qty	Case Qty
73230-00	5cm x 7cm	1 x 10	28 x 10
73230-01	7cm x 10cm	1 x 10	20 x 10
73230-02	8cm x 15cm	1 x 10	20 x 10
73230-03	10cm x 24cm	1 x 5	20 x 5
73230-05	10cm x 35cm	1 x 5	10 x 5



Leukomed®

Non-woven dressing

- A non-woven dressing, gently adherent with integrated absorbent wound pad
- Highly breathable and conformable dressing

Ideal for: surgical or traumatic wounds with low to moderate exudate levels



Code	Size	Pack Qty	Case Qty
72380-00	7.2cm x 5cm	1 x 50	8 x 50
72380-01	8cm x 10cm	1 x 50	6 x 50
72380-02	8cm x 15cm	1 x 50	6 x 50
72380-03	10cm x 20cm	1 x 50	8 x 50
72380-04	10cm x 25cm	1 x 50	8 x 50
72380-05	10cm x 30cm	1 x 50	8 x 50
72380-06	10cm x 35cm	1 x 50	8 x 50

Leukomed® T

Transparent film dressing

- Ultra thin and conformable transparent film dressing
- Forms a sterile barrier reducing the risk of wound infection
- Allows easy wound inspection without need to remove dressing

Ideal for: sterile fixation of dressings and instruments and protection of dry wounds



Code	Size	Pack Qty	Case Qty
72381-00	7.2cm x 5cm	1 x 50	10 x 50
72381-01	8cm x 10cm	1 x 50	6 x 50
72381-08	10cm x 12.5cm	1 x 50	8 x 50
72381-14	10cm x 25cm	1 x 50	8 x 50
72381-02	11cm x 14cm	1 x 50	8 x 50
72381-09	15cm x 20cm	1 x 50	8 x 50
72381-10	15cm x 25cm	1 x 50	8 x 50

Leukomed® T plus

Transparent film dressing with wound pad

- Ultra thin and conformable film dressing with a low-adherent wound pad
- Sterile barrier protects the wound and supports undisturbed wound healing

Ideal for: dressing surgical incisions, cover and protect superficial cuts, lacerations and minor burns



Code	Size	Pack Qty	Case Qty
72382-00	7.2cm x 5cm	1 x 50	8 x 50
72382-01	8cm x 10cm	1 x 50	6 x 50
72382-02	8cm x 15cm	1 x 50	6 x 50
72382-14	10cm x 20cm	1 x 50	8 x 50
72382-03	10cm x 25cm	1 x 50	8 x 50
72382-04	10cm x 30cm	1 x 50	8 x 50
72382-05	10cm x 35cm	1 x 50	8 x 50

Cuticell® range

Cuticell® Classic

Paraffin-impregnated gauze dressing

- A low-adherent dressing made from open weave cotton gauze, impregnated with soft paraffin
- Supports a moist wound bed environment

Ideal for: Exuding wounds, including first and second degree burns, lacerations, leg ulcers, donor and recipient skin graft sites

Code	Size	Pack Qty	Case Qty
Cuticell Classic			
72538-00	5cm x 5cm - Individually sealed	1 x 5	120 x 5
72538-01	5cm x 5cm - Individually sealed	1 x 50	12 x 50
72538-02	10cm x 10cm - Individually sealed	1 x 10	60 x 10
72538-03	10cm x 10cm - Individually sealed	1 x 100	6 x 100
72538-04	10cm x 40cm - Individually sealed	1 x 10	30 x 10
72538-05	10cm x 10cm - Plastic box	1 x 36	12 x 36
72538-06	10cm x 7m - Dispenser roll	1 x 1	1 x 12

Code	Size	Pack Qty	Case Qty
Cuticell Classic			
72680-00	5cm x 7.5cm	1 x 5	-
72680-01	7.5cm x 10cm	1 x 5	-
72680-02	10cm x 18cm	1 x 5	-
72680-03	15cm x 25cm	1 x 5	-

Cuticell® Contact

Silicone wound contact layer

- Gentle and highly conformable silicone wound contact layer for pain-free and atraumatic dressing changes



Soffban® Synthetic

Synthetic padding

- Made from 100% polyester
- Autoclavable
- Soft and highly conformable
- Maximum loft retention for excellent cushioning
- Reduced risk of skin maceration and sensitivity
- Tearability for easy application

Ideal for: use under synthetic casts, plaster of Paris and bandages



Code	Size	Pack Qty	Case Qty
71464-01	5cm x 2.7m	1 x 12	24 x 12
71467-01	7.5cm x 2.7m	1 x 12	16 x 12
71486-05	10cm x 2.7m	1 x 12	12 x 12
71465-01	15cm x 2.7m	1 x 12	8 x 12
71466-01	20cm x 2.7m	1 x 6	12 x 6

Soffban® Natural

Natural padding

- Suitable to use where bleeding can occur
- A highly absorbent padding made from 100% viscose rayon
- Excellent conformability
- Easy to tear
- Autoclavable

Ideal for: using as a comfortable undercast layer, use under bandages or under plaster of Paris and synthetic casts



Code	Size	Pack Qty	Case Qty
71472-00	5cm x 2.7m	1 x 12	24 x 12
71473-00	7.5cm x 2.7m	1 x 12	16 x 12
71469-00	10cm x 2.7m	1 x 12	12 x 12
71470-00	15cm x 2.7m	1 x 12	8 x 12
71471-00	20cm x 2.7m	1 x 6	12 x 6

Easifix®

Retention bandage

- Made from a quality viscose and polyamide material to provide a soft, lightweight and porous bandage
- Excellent conformability

Ideal for: fixation of dressings on all body parts especially joints



Code	Size	Pack Qty	Case Qty
71429-10	5cm x 4m stretched	1 x 20	24 x 20
71430-03	7.5cm x 4m stretched	1 x 20	24 x 20
71428-04	10cm x 4m stretched	1 x 20	12 x 20
71429-11	15cm x 4m stretched	1 x 20	12 x 20

Easifix® K

Knitted retention bandage

- Knitted bandage made from viscose and polyamide, giving two-way stretch
- Lightweight, soft and breathable
- Highly elastic and flexible bandage that easily conforms to the most awkward body parts

Ideal for: fixation of dressings on all body parts especially joints



Code	Size	Pack Qty	Case Qty
72617-00	2.5cm x 4m	1 x 20	24 x 20
72617-01	5cm x 4m	1 x 20	24 x 20
72617-02	7.5cm x 4m	1 x 20	24 x 20
72617-03	10cm x 4m	1 x 20	12 x 20
72617-04	15cm x 4m	1 x 20	12 x 20



Tensoplast®

Elastic adhesive bandage


- High quality bandage made from a combination of specially formulated porous adhesive and a superior woven cloth
- Soft, non-fraying edges
- Standard and reverse wound bandage option
- Reverse wound bandage has adhesive part facing outwards



Ideal for: quick, easy and safe application to horses hind legs

Code	Size	Pack Qty	Case Qty
72050-58	2.5cm x 4.5m (plaster)	1 x 12	12 x 12
72050-32	5cm x 4.5m stretched	1 x 12	6 x 12
72050-36	7.5cm x 4.5m stretched	1 x 12	6 x 12
72050-29	10cm x 4.5m stretched	1 x 12	6 x 12
72050-59	5cm x 4.5m stretched (reverse wound)	1 x 60	1 x 60
72050-34	7.5cm x 4.5m stretched (reverse wound)	1 x 48	1 x 48

Tensoplast® Reverse Wound Elastic adhesive bandage



Tensoplast Reverse Wound bandage is a preferred EAB when bandaging the hind legs of a difficult horse. The adhesive part facing outwards allows for easy roll on application and with little drag effect, making it quicker and therefore safer to apply. This also gives more control with regards to how much tension is applied.



Easifix® Cohesive

Cohesive retention bandage

- Latex-coated cohesive bandage made from viscose and polyester
- Securely sticks to itself for strong fixation of wound dressings and plaster splints
- Only a 30% overlap is required



Ideal for: secure fixation of wound dressings

Code	Size	Pack Qty	Case Qty
71436-14	4cm x 20m stretched	1 x 1	66 x 1
71436-15	6cm x 20m stretched	1 x 1	1 x 48
71436-16	8cm x 20m stretched	1 x 1	1 x 36
71436-17	10cm x 20m stretched	1 x 1	1 x 30
71436-18	12cm x 20m stretched	1 x 1	1 x 24
71436-09	4cm x 4m stretched	1 x 10	300 x 10
71436-10	6cm x 4m stretched	1 x 10	240 x 10
71436-11	8cm x 4m stretched	1 x 10	180 x 10
71436-12	10cm x 4m stretched	1 x 10	150 x 10
71436-13	12cm x 4m stretched	1 x 10	120 x 10

Co-Plus®

Cohesive bandage

- A non-slip, light support bandage suitable for final layer fixation
- Flexible bandage, conforms easily to the most awkward body parts
- Holds in place by sticking to itself

Ideal for: light compression or as secondary fixation



Code	Size	Pack Qty	Case Qty
White			
72101-09	5cm x 6.3m stretched	1 x 36	1 x 36
72101-10	7.5cm x 6.3m stretched	1 x 24	1 x 24
72101-05	10cm x 6.3m stretched	1 x 18	1 x 18
Mixed			
72102-04	5cm x 6.3m stretched	1 x 36	1 x 36
72102-03	7.5cm x 6.3m stretched	1 x 24	1 x 24
72102-02	10cm x 6.3m stretched	1 x 18	1 x 18
Dark blue			
72102-23	10cm x 6.3m stretched	1 x 18	1 x 18
Green			
72102-05	10cm x 6.3m stretched	1 x 18	1 x 18

Delta Terry-Net® C

Natrual Terry Cloth Stockinette

- Cotton, soft stockinette, which provides good cushioning and protection
- Easy to use dispenser
- Autoclavable

Ideal for: underneath padding or synthetic and plaster of Paris casting materials



Code	Size	Pack Qty	Case Qty
72305-00	5cm x 13.7m	1 x 1	1 x 2
72305-01	7.5cm x 13.7m	1 x 1	1 x 2
72305-02	10cm x 13.7m	1 x 1	1 x 2
72305-03	12.5cm x 13.7m	1 x 1	1 x 1
72305-04	15cm x 13.7m	1 x 1	1 x 1



Delta-Cast® Elite

Non-fibreglass cast tape

- An extremely strong and lightweight cast tape made from polyester substrate
- Extensible yarns which shape easily to body contours
- Soft edges and smooth finish



Ideal for: primary and secondary cast requirements for both small and large animals

Code	Size	Pack Qty	Case Qty
72279-00	5cm x 3.6m	1 x 10	1 x 10
72279-01	7.5cm x 3.6m	1 x 10	1 x 10
72279-02	10cm x 3.6m	1 x 10	1 x 10
72279-03	12.5cm x 3.6m	1 x 10	1 x 10

Delta-Cast® Conformable

Non-Fibreglass Cast Tape

Delta-Cast® Conformable has set a benchmark in standard and removable, non-fibreglass cast tapes. Delta-Cast® Conformable responds to the daily challenges of fracture management - it's distinctive construction provides rigid and semi-rigid casting options for primary and secondary casting applications where patient fit, comfort and compliance are critical.

Delta-Cast® Conformable is the only UK cast tape with clinical evidence¹ to support its use for focused rigidity casting.



Code	Size	Colour	Pack Qty
72280-00	2.5cm x 1.80m	○ White	1 x 10
72280-01	5cm x 3.6m	○ White	1 x 10
72280-02	7.5cm x 3.6m	○ White	1 x 10
72280-03	10cm x 3.6m	○ White	1 x 10
72280-04	12.5cm x 3.6m	○ White	1 x 10

1. Petty, A. C., Wardman, C., (1998). A randomized, controlled comparison of adjustable focused rigidity primary casting technique with standard plaster of Paris/ synthetic casting technique in the management of fractures and other injuries. Journal of Orthopaedic Nursing, 2, 95 - 102

Gypsona®

Premium plaster of Paris

- Plaster of Paris bandage providing excellent stability and conformability using an interlocked leno weave gauze
- High quality, fast setting plaster of Paris
- Available as a roll or slab



Ideal for: supporting and immobilising fresh fractures and soft tissue injuries

Code	Size	Pack Qty	Case Qty
71988-26	90cm x 5.5m roll	1 x 1	1 x 1
71985-04	5cm x 2.7m roll	1 x 48	1 x 48
71985-05	7.5cm x 2.7m roll	1 x 48	1 x 48
71985-00	10cm x 2.7m roll	1 x 48	1 x 48
71985-02	15cm x 2.7m roll	1 x 24	1 x 24
71985-03	20cm x 2.7m roll	1 x 24	1 x 24
71988-14	5 ply slab 10cm x 20m	1 x 1	1 x 1
71988-15	5 ply slab 15cm x 20m	1 x 1	1 x 1
71988-16	5 ply slab 20cm x 20m	1 x 1	1 x 1

Dynacast® Prelude

Synthetic splint system

- A strong, light and durable synthetic splint system incorporating fibreglass, covered with a polypropylene padding
- Quick setting and allows early weight bearing



Ideal for: supporting and immobilising fresh fractures and soft tissue injuries

Code	Size	Pack Qty	Case Qty
71418-00	2.5cm x 4.6m	1 x 1	2 x 1
71420-00	5cm x 4.6m	1 x 1	2 x 1
71421-00	7.5cm x 4.6m	1 x 1	2 x 1
71415-00	10cm x 4.6m	1 x 1	2 x 1
71416-00	12.5cm x 4.6m	1 x 1	2 x 1
71417-00	15cm x 4.6m	1 x 1	2 x 1
71419-00	20cm x 2.3m	1 x 1	2 x 1

Hypafix®

Wide-area fixation

- Self-adhesive, non-woven wide area dressing fixation
- Easily cut to shape and size required

Ideal for: initial and long term primary dressing adhesion



Code	Size	Pack Qty	Case Qty
71442-03	5cm x 5m	1 x 1	24 rolls
71442-04	10cm x 5m	1 x 1	12 rolls
71443-00	2.5cm x 10m	1 x 24	24 rolls
71443-01	5cm x 10m	1 x 1	24 rolls
71443-02	10cm x 10m	1 x 1	24 rolls
71443-03	15cm x 10m	1 x 1	24 rolls
71443-04	20cm x 10m	1 x 1	24 rolls
71443-05	30cm x 10m	1 x 1	24 rolls

Leukopor®

Gentle tape

- A gentle fixation tape made from a non-woven polyester viscose fabric
- Good initial and long-term adhesion
- Can be torn easily by hand

Ideal for: sensitive skin and fixation of dressings and instruments



Code	Size	Pack Qty	Case Qty
02471-00	1.25cm x 5m	1 x 24	10 x 24
02472-00	2.5cm x 5m	1 x 12	10 x 12
02474-00	5cm x 5m	1 x 6	10 x 6
02453-00	1.25cm x 9.2m	1 x 24	10 x 24
02454-00	2.5cm x 9.2m	1 x 12	10 x 12
02455-00	5cm x 9.2m	1 x 6	10 x 6

Hypafix® Transparent

Waterproof dressing fixation

- Thin and conformable film
- Protects dressings from contamination
- Excellent adhesion and flexibility

Ideal for: waterproof wide-area dressing fixation



Code	Size	Pack Qty	Case Qty
72378-00	5cm x 10m	1 x 1	1 x 12
72378-01	10cm x 10cm	1 x 1	1 x 12
72378-02	15cm x 10m	1 x 1	1 x 12
72378-03	10cm x 2m	1 x 1	1 x 12

Leukosilk®

The ideal tape for sensitive skin

- The high strength universal tape made from a fine acetate fabric
- Good initial and permanent adhesion
- Easy to write on

Ideal for: applications such as securing tubes and other instruments



Code	Size	Pack Qty	Case Qty
72669-00	1.25cm x 4.6m	1 x 24	10 x 24
72669-01	2.5cm x 4.6m	1 x 12	10 x 12
72669-02	5cm x 4.6m	1 x 6	10 x 6

Leukofix®

Transparent tape

- Skin-friendly transparent polyethylene film with polyacrylate adhesive
- Provides secure fixation and visibility of the puncture site
- Fits all body contours

Ideal for: fixation of dressings and instruments of all types and sizes



Code	Size	Pack Qty	Case Qty
02121-00	1.25cm x 5m	1 x 24	10 x 24
02122-00	2.5cm x 5m	1 x 12	10 x 12
02124-00	5cm x 5m	1 x 6	10 x 6
02136-00	1.25cm x 9.2m	1 x 24	10 x 24
02137-00	2.5cm x 9.2m	1 x 12	10 x 12
02138-00	5cm x 9.2m	1 x 6	10 x 6

Leukoplast® Sleek

Strong, waterproof adhesive tapes

- A smooth, thin, plastic film
- Strong, waterproof and reliable
- Latex-free variants available

Ideal for: securely holding dressings and instruments in place



Code	Size	Pack Qty	Case Qty
72379-01	2.5cm x 3m	1 x 12	18 x 12
72379-02	2.5cm x 5m	1 x 12	18 x 12
72379-03	5cm x 5m	1 x 12	12 x 12
72379-04	7.5cm x 5m	1 x 12	8 x 12
72379-05	10.5cm x 5m	1 x 24	1 x 24

Tensoplast®

Elastic adhesive bandage

- High quality bandage made from a combination of specially formulated porous adhesive and a superior woven cloth
- Soft, non-fraying edges
- Standard and reverse wound bandage option
- Reverse wound bandage has adhesive part facing outwards

Ideal for: quick, easy and safe application to horses hind legs



Code	Size	Pack Qty	Case Qty
72050-58	2.5cm x 4.5m (plaster)	1 x 12	12 x 12

Strappal®

Rigid strapping tape

- Constructed from a viscose fabric and zinc oxide adhesive
- Hand tearable for quick and secure application

Ideal for: strong and effective stabilisation of joint structures



Code	Size	Pack Qty	Case Qty
71496-15	1.25cm x 10m	1 x 36	1 x 36
71487-00	2.5cm x 10m	1 x 36	1 x 36
71490-00	4cm x 10m	1 x 24	1 x 24
71493-00	5cm x 10m	1 x 18	1 x 18
71494-01	7.5cm x 10m	1 x 6	1 x 6

Propax® Gauze Swabs

Gauze swabs

- Made from high quality BP Type 13 gauze
- Available in two colours; green for anaesthetics and white for general use
- Available sterile and non-sterile

Ideal for: wound cleansing, skin preparation, general cleansing or as a primary or secondary dressing

Code	Size	Pack Qty	Case Qty
Propax Gauze swabs Type 13 BP (Sterile)			
71117-02	7.5cm x 7.5cm 5 pads - 8 ply	1 x 48	12 x 48
Propax Gauze swabs Type 13 BP (Non-sterile)			
71734-06	5cm x 5cm - 8 ply	1 x 100	50 x 100
71734-05	5cm x 5cm - 12 ply	1 x 100	50 x 100
71734-12	7.5cm x 7.5cm - 8 ply	1 x 100	50 x 100
71734-13	7.5cm x 7.5cm - 12 ply	1 x 100	48 x 100
71734-10	10cm x 10cm - 8 ply	1 x 100	50 x 100
71734-11	10cm x 10cm - 12 ply	1 x 100	50 x 100
Propax Gauze swabs Type 13 BP (Non-sterile) - Green			
71734-08	7.5cm x 7.5cm - 8 ply	50 x 100	50 x 100
71734-07	10cm x 10cm - 8 ply	1 x 100	50 x 100
71736-28	10cm x 10cm - 12 ply	40 x 100	40 x 100



Multisorb®

Non-woven swabs

- Low linting, soft and absorbent swabs
- Available sterile or non-sterile

Ideal for: cleaning wounds, skin preparation or as a secondary dressing

Code	Size	Pack Qty	Case Qty
Multisorb non-woven swabs (Sterile)			
71786-05	5cm x 5cm - 4 ply	50 x 5	12 x 50 x 5
71787-08	7.5cm x 7.5cm - 4 ply	50 x 5	12 x 50 x 5
71786-07	10cm x 10cm - 4 ply	50 x 5	8 x 50 x 5
Multisorb non-woven swabs (Non-sterile)			
71783-06	7.5cm x 7.5cm - 6 ply	50 x 50	50 x 50
71783-04	10cm x 10cm - 6 ply	50 x 50	50 x 50



Fivex®

X-Ray detectable swabs

- Made from high quality, BP Type 13 gauze
- Available in two colours; green for anaesthetics and white for general use
- Available sterile and non-sterile

Ideal for: surgical procedures

Code	Size	Pack Qty	Case Qty
71766-20	15cm x 2.5cm - 16 ply	10 x 10 T5s	10 x 10 T5s

A Quick Summary Guide to our Product Range

Wound Dressings	Bandages	Tapes	Casting
<p>Wound bed management</p> <p>Cutimed Gel</p> <ul style="list-style-type: none"> Amorphous hydrogel for necrotic or sloughy wounds or skin hydration <p>Cutimed Debriclean</p> <ul style="list-style-type: none"> Mechanical Debridement Pad 	<p>Elastic adhesive bandage</p> <p>Tensoplast</p> <ul style="list-style-type: none"> High quality bandage made from a combination of specially formulated porous adhesive and a superior woven cloth <p>Tensoplast (Reverse)</p> <ul style="list-style-type: none"> Easy roll on for quicker and therefore safer application. Preferred EAB for bandaging hind legs of a difficult horse 	<p>Low-allergy tape</p> <p>Leukopor</p> <ul style="list-style-type: none"> Suitable for long term dressing fixation on sensitive skin 	<p>Casting materials</p> <p>Delta-Cast Elite</p> <ul style="list-style-type: none"> An extremely strong and lightweight cast tape made from polyester substrate <p>Delta-Cast Conformable</p> <ul style="list-style-type: none"> Removable, non-fibreglass cast tapes. Providing rigid and semi-rigid casting options for primary and secondary casting applications where patient fit, comfort and compliance are critical <p>Gypsona</p> <ul style="list-style-type: none"> Plaster of Paris bandage providing excellent stability and conformability using an interlocked leno weave gauze <p>Dynacast Prelude</p> <ul style="list-style-type: none"> A strong, light and durable synthetic splint system incorporating fibreglass, covered with a polypropylene padding
<p>Exudate management</p> <p>Cutimed Siltec</p> <ul style="list-style-type: none"> Foam dressings with a silicone wound contact layer and super-absorbent technology for effective fluid management <p>Cutimed Sorbion Sachet Extra</p> <ul style="list-style-type: none"> Super-absorbent dressings for high levels of exudate 	<p>Cohesive bandage</p> <p>Co-Plus</p> <ul style="list-style-type: none"> Adheres only to itself. Fixation as final layer 	<p>Waterproof adhesive tape</p> <p>Leukoplast Sleek</p> <ul style="list-style-type: none"> For secure fixation where contact with dirt or water is likely 	<p>Swabs</p> <p>Swabs</p> <p>Propax</p> <ul style="list-style-type: none"> Gauze swabs made from high quality BP Type 13 <p>Multisorb</p> <ul style="list-style-type: none"> Non-woven, low linting, soft and absorbent swabs <p>Fivex</p> <ul style="list-style-type: none"> X-Ray detectable swabs made from high quality BP Type 13
<p>Infection management</p> <p>Cutimed Sorbact</p> <ul style="list-style-type: none"> Irreversibly binds bacteria to its surface to reduce bioburden and support wound healing 	<p>Synthetic padding</p> <p>Soffban Synthetic</p> <ul style="list-style-type: none"> Providing cushioning and protection of injuries 	<p>Transparent tape</p> <p>Leukofix</p> <ul style="list-style-type: none"> Fixation of dressings, canulae, probes and catheters 	
<p>Tissue management</p> <p>Cutimed Epiona</p> <ul style="list-style-type: none"> Boosts healing in chronic and stagnating wounds 	<p>Knitted retention bandage</p> <p>Easifix K</p> <ul style="list-style-type: none"> Highly conformable and breathable open-knitted bandage with two-way stretch 	<p>Wide area fixation</p> <p>Hypafix</p> <ul style="list-style-type: none"> Suitable for the fixation of dressings Easily cut to the correct size and shape 	
<p>Post-operative</p> <p>Leukomed Sorbact</p> <ul style="list-style-type: none"> Innovative post-operative dressing to help to reduce wound bacteria colonisation by its physical mode of action. For post-operative and traumatic wounds 	<p>Cohesive retention bandage</p> <p>Easifix Cohesive</p> <ul style="list-style-type: none"> For secure fixation of wound dressings 	<p>Rigid strapping tape</p> <p>Strappal</p> <ul style="list-style-type: none"> Strong support and efficient stabilisation of joint structures. Hypoallergenic 	
<p>Other</p> <p>Leukomed</p> <ul style="list-style-type: none"> Non-woven island dressing with an absorbent wound pad Suitable for post-op and minor injuries <p>Leukomed T plus</p> <ul style="list-style-type: none"> Waterproof dressing suitable for post-op and minor injuries 			

Education by Essity

The Essity veterinary team are excited to be able to conduct wound care education and product training with you via online meetings.

Training topics covered include;

- General wound care products
- Debridement
- Exudate management
- Post-operative wound protection

We can tailor the training to suit your needs.

Book a training session with us and you will receive a free advanced wound care sample pack which will contain items specific to your training.

Contact **Francesca Gordon, animalhealthcare@essity.com**

Please contact the team for more information

Tel: **01482 670124**

Email: **animalhealthcare@essity.com**