

Cutimed® PROTECT

Think Skin
Protection.
Think Cutimed®
PROTECT



Cutimed PROTECT cream - 2g sachet now available

Cutimed®

Cutimed PROTECT

Medical skin protection

The skin is the human body's largest organ. Its basic function is to provide protection from external factors as well as prevention of trans-epidermal water loss by maintaining a protective barrier. Cutimed PROTECT can help to re-establish and maintain the skin's barrier function.



Cutimed PROTECT film spray and foam applicator

Supplied as a liquid Cutimed PROTECT film forms a transparent, durable and waterproof barrier that works as a protective shield between the skin and bodily waste fluids, adhesive products and shearing forces. Cutimed PROTECT spray and foam applicator offer a long-lasting protective barrier for up to 96 hours.

- Protects wound margins and peri-wound skin from external factors including wound exudate, incontinence and friction
- Preserves and maintains intact or irritated skin by building a moisture barrier with good permeability to oxygen and water vapour without loss of adhesion for dressings
- Prevents maceration by providing up to 96 hours of moisture barrier
- Quick drying, breathable and alcohol-free the liquid does not sting on application and can be removed by mild medical adhesive remover if necessary

Cutimed PROTECT cream

As a barrier cream Cutimed PROTECT preserves and maintains intact skin by providing a protective barrier that can allow adhesive products, such as dressings, to adhere to the skin as normal.

- Moisturises the skin and prevents skin from drying out while keeping it permeable to water vapour and oxygen
- Preventive treatment against maceration caused by incontinence or body fluids
- Quickly and easily rubbed in with little friction and without leaving a greasy film, it is alcohol- and perfume-free

Code	Size	Pack Qty	PIP Code	NHS SC Code
Cream - For protection of skin at risk from damage and irritation				
72652-08	2g	1 x 20	405-1447	ELY614
72652-00	28g	1 x 1	369-4916	ELY471
72652-03	90g	1 x 12	369-4924	ELY472
Spray - For wide area application				
72653-00	28ml	1 x 1	369-4882	-
72652-08	28ml	1 x 12	-	-
Foam Applicator - For easy and precise application around wound margins				
72654-00	1ml	1 x 5	369-4890	-
72654-01	1ml	1 x 25	-	-
72654-02	3ml	1 x 5	369-4908	-
72654-03	3ml	1 x 25	-	-

Contraindications:

Cutimed PROTECT cream should not be used on skin where infection is present or when known sensitivity to the ingredients.

To arrange a visit from your local Account Manager or order a free evaluation sample, please contact the **Concierge Service** on **01482 670177**