

# JOBST® UltraSheer


Ready-to-Wear RAL Compression Hosiery

Available on Drug Tariff




Comfort, Health and Style!

RAL Compression Classes 1 (18-21mmHg) and 2 (23-32mmHg)						
Measuring point	I	II	III	IV	V	VI
cG silicone band	43-48	45-52	49-56	53-60	56-64	60-68
cG silicone band wide	49-57	53-62	57-67	61-72	65-77	69-81
cG	43-57	45-62	49-67	53-72	56-77	60-81
cF	39-52	41-56	44-60	47-65	50-69	53-73
cE	30-37	33-40	35-43	37-45	39-48	41-51
cD	27-33	29-36	32-39	34-42	36-45	38-48
cC	28-34	30-37	33-40	35-43	37-46	39-49
cB1	23-27	24-29	26-32	29-35	31-37	33-39
cB	18-20	20-22	22-24	24-26	26-28	28-30
cY	26-31	28-33	29-35	31-37	32-38	33-40
cA	17-22	19-24	21-26	23-29	25-32	27-34
Length measurements for all sizes in cm						
AD: ID regular	39-45cm			ID petite	34-39cm	
AG: IG regular	70-80cm			IG petite	60-69cm	
AT: IK regular	75-85cm			IK petite	65-75cm	



Please note: not all colours and options available for all styles and compression classes.

Customer Services: **0845 122 3600**



### Silicone Top Band Thigh High Options

