

Available on prescription

Custom-fit, flat-knit compression garments

STRONG IN ACTION AND PERFORMANCE


Knee High (AD)

Basic Style - slant open toe in beige

- L1-01-04 CCL1 - 18-21mmHg
- L2-02-04 CCL2 - 23-32mmHg
- L3-03-04 CCL3 - 34-46mmHg
- L3-04-04 CCL3 Forte - 34-46mmHg
- L4-05-04 CCL4 - 49-70mmHg
- L5-06-04 CCL4 Super - 60-90mmHg

Options

- L-A001 Closed toe
- L-A002 2 Ankle pads (profile)
- L-A003 Zipper
- L-A004 Silicone band
- L-A008 Colour option: black
- L-A010 T-heel



Tights (AT)

Basic Style - slant open toe in beige

- L1-01-06 CCL1 - 18-21mmHg
- L2-02-06 CCL2 - 23-32mmHg
- L3-03-06 CCL3 - 34-46mmHg
- L3-04-06 CCL3 Forte - 34-46mmHg
- L4-05-06 CCL4 - 49-70mmHg
- L5-06-06 CCL4 Super - 60-90mmHg

Options

- L-A001 Closed toe
- L-A002 2 Ankle pads (profile)
- L-A003 Zipper
- L-A004 Silicone band
- L-A007 Leg extension
- L-A008 Colour option: black
- L-A009 Slip-form
- L-A0010 T-heel




Gauntlet to Elbow (AE)

Basic Style - no fingers, no thumb in beige

- U1-01-03 CCL1 - 14-18mmHg
- U2-02-03 CCL2 - 20-25mmHg
- U3-03-03 CCL3 - 25-30mmHg

Options

- U-A001 Zipper
- U-A002 Silicone band
- U-A007 Colour option: black
- U1-A009 Finger - CCL1 - 14-18mmHg
- U2-A011 Finger - CCL2 - 20-25mmHg
- U3-A013 Finger - CCL3 - 25-30mmHg




Thigh High (AG)

Basic Style - slant open toe in beige

- L1-01-06 CCL1 - 18-21mmHg
- L2-02-06 CCL2 - 23-32mmHg
- L3-03-06 CCL3 - 34-46mmHg
- L3-04-06 CCL3 Forte - 34-46mmHg
- L4-05-06 CCL4 - 49-70mmHg
- L5-06-06 CCL4 Super - 60-90mmHg

Options

- L-A001 Closed toe
- L-A002 2 Ankle pads (profile)
- L-A003 Zipper
- L-A004 Silicone band
- L-A008 Colour option: black
- L-A009 Slip-form
- L-A010 T-heel



Foot Cap with Toes (AA)

Beige

- L1-01-02 CCL1 - 18-21mmHg
- L2-02-02 CCL2 - 23-32mmHg
- L3-03-02 CCL3 - 34-46mmHg

Options

- L-A008 Colour option: black



Armsleeve to Axilla with Gauntlet (AG)

Basic Style - no fingers, no thumb in beige

- U1-01-05 CCL1 - 14-18mmHg
- U2-02-05 CCL2 - 20-25mmHg
- U3-03-05 CCL3 - 25-30mmHg

Options

- U-A001 Zipper
- U-A002 Silicone band
- U-A003 Shoulder cap
- U-A007 Colour option: black
- U1-A009 Finger - CCL1 - 14-18mmHg
- U2-A011 Finger - CCL2 - 20-25mmHg
- U3-A013 Finger - CCL3 - 25-30mmHg




Chap Style (AG-T) Thigh High with Waist Attachment

Basic Style - slant open toe in beige

- L1-01-06 CCL1 - 18-21mmHg
- L2-02-06 CCL2 - 23-32mmHg
- L3-03-06 CCL3 - 34-46mmHg
- L3-04-06 CCL3 Forte - 34-46mmHg
- L4-05-06 CCL4 - 49-70mmHg
- L5-06-06 CCL4 Super - 60-90mmHg

Options

- L-A001 Closed toe
- L-A002 2 Ankle pads (profile)
- L-A003 Zipper
- L-A004 Silicone band
- L-A008 Colour option: black
- L-A009 Slip-form
- L-A010 T-heel



Gauntlet to Wrist (AC)

Basic Style - no fingers, no thumb in beige

- U1-01-01 CCL1 - 14-18mmHg
- U2-02-01 CCL2 - 20-25mmHg
- U3-03-01 CCL3 - 25-30mmHg

Options

- U-A001 Zipper
- U-A002 Silicone band
- U-A007 Colour option: black
- U1-A009 Finger - CCL1 - 14-18mmHg
- U1-A011 Finger - CCL2 - 20-25mmHg
- U1-A013 Finger - CCL3 - 25-30mmHg




Armsleeve from Wrist to Axilla (CG)

Basic Style - in beige

- U1-01-07 CCL1 - 14-18mmHg
- U2-02-07 CCL2 - 20-25mmHg
- U3-03-07 CCL3 - 25-30mmHg

Options

- U-A001 Zipper
- U-A002 Silicone band
- U-A003 Shoulder cap
- U-A007 Colour option: black



Lymphoedema Helpline:

Freephone 0800 389 8424

or email:

lymphoedema.uk@bsnmedical.com



Compression Classes

Lower limb garments (RAL)

CCL1	CCL2	CCL3	CCL3 Forte	CCL4	CCL4 Super
18-21 mmHg	23-32 mmHg	34-46 mmHg	34-46 mmHg	49-70 mmHg	60-90 mmHg

Upper limb garments

CCL1	CCL2	CCL3
14-18 mmHg	20-25 mmHg	25-30 mmHg

Notes: CCL3 Forte, 4 and 4 Super are available for the lower body only.

Special styles are available for most custom-fit garments e.g. open pubis, zips, adjustable waist bands and fly for men (most available on Drug Tariff).