

JOBST® Opaque



Soft and comfortable leg therapy – all day long!

Measurement Point	RAL Compression classes 1 (18-21mmHg) and 2 (23-32mmHg)					
	I	II	III	IV	V	VI
cG silicone band	43-48	45-52	49-56	53-60	56-64	60-68
cG silicone band wide	49-57	53-62	57-67	61-72	65-77	69-81
cG	43-57	45-62	49-67	53-72	56-77	60-81
cF	39-52	41-56	44-60	47-65	50-69	53-73
cE	30-37	33-40	35-43	37-45	39-48	41-51
cD	27-33	29-36	32-39	34-42	36-45	38-48
cC	28-34	30-37	33-40	35-43	37-46	39-49
cB1	23-27	24-29	26-32	29-35	31-37	33-39
cB	18-20	20-22	22-24	24-26	26-28	28-30
cY	26-31	28-33	29-35	31-37	32-38	33-40
cA	17-22	19-24	21-26	23-29	25-32	27-34

Length measurements for all sizes (in cm):

AD: ID long 39-45, ID short 34-39 AG: IG long 70-80, IG short 60-69 AT/AG-T: IK long 75-85, IK short 65-75

Ready-to-Wear
RAL Compression
Hosiery



Circular-knit, seamless
compression hosiery
manufactured using
JOBST® Advanced Comfort

- Softer than ever for exceptional wearing comfort
- Effective moisture wicking for well tempered legs
- Great fit with superior donning

RAL Compression

- Compression class 1 (18-21mmHg)
- Compression class 2 (23-32mmHg)


Five Colours



Six Size Options

SIZE					
I	II	III	IV	V	VI
LONG					
SHORT					


Garment Styles



AD Knee High*

AG Thigh High*


AG-T Chap




AT Tights*

AT Men

AT Maternity



Top Band Thigh High Options



Knitted top band

Soft silicone top band*

Dotted silicone top band*

Lace silicone top band*

*These styles/colours available on prescription



BSN medical Limited.
PO Box 258, Willerby, Hull HU10 6WT
Tel: 01482 670100 Fax: 01482 670111
www.bsnmedical.co.uk

www.jobst.com
JOBST - a brand of **BSN medical**

For product and ordering information call **0845 122 3600** (Option 1)
Customer Services Direct Fax: **0845 122 3450**
Lymphoedema Helpline: **0800 389 8424**
Email: lymphoedema.uk@bsnmedical.com

JOBST®
Comfort, Health and Style!