

JOBST® Confidence

Custom-Fit Order Form

Fax order to customer services on: **0845 122 3450**

Email order to customer services on: **compression.uk@jobst.com**

Date: _____ Purchase Order No.: _____ Patient Name: _____ DoB: _____
 Measured By: _____ Tel: _____ Email: _____
 Delivery Address: _____ Invoice Address: _____



<input type="checkbox"/> JOBST® Confidence	
Compression Class (RAL)	Quantity
<input type="checkbox"/> CCL 1 (18-21mmHg)	Left <input type="text"/>
<input type="checkbox"/> CCL 2 (23-32mmHg)	Right <input type="text"/>
<input type="checkbox"/> CCL 3 (34-46mmHg)	
Style and Options	
<input type="checkbox"/> AD Knee high	
<input type="checkbox"/> Slipform	<input type="checkbox"/> Ankle functional zone
<input type="checkbox"/> Decor-line	<input type="checkbox"/> Initials (max 2 letters) <input type="text"/>
<input type="checkbox"/> Hallux valgus (bunion relief zone)*	
<input type="checkbox"/> AG Thigh high	
<input checked="" type="checkbox"/> Slipform	<input type="checkbox"/> Ankle functional zone
<input type="checkbox"/> Decor-line	<input type="checkbox"/> Knee functional zone
<input type="checkbox"/> Initials (max 2 letters) <input type="text"/>	<input type="checkbox"/> Hallux valgus (bunion relief zone)*
Colour	Top band options
<input type="checkbox"/> Beige	<input type="checkbox"/> Black
<input type="checkbox"/> Caramel	<input type="checkbox"/> Red heather
<input type="checkbox"/> Anthracite heather	<input type="checkbox"/> Jeans heather
Note: the functional zone for red, anthracite and jeans heather coloured garments is more visible	
Knee high A-D	
<input type="checkbox"/> 5cm SoftFit	
Thigh high A-G	
<input checked="" type="checkbox"/> 5cm silicone top band	
Remarks	

Left Leg Length Measurements (cms)

a-K _____
 Required for thigh high garment

a-G _____ G _____

a-F _____ F _____

a-E _____ E _____

a-D _____ **D** _____
 (Knee high)

a-C _____ C _____


a-B' _____ B' _____

a-B _____ B _____

_____ Y _____

_____ A _____

Circumference Measurements (cms)



Right Leg Length Measurements (cms)

a-K _____
 Required for thigh high garment

G _____ a-G _____

F _____ a-F _____

E _____ a-E _____

D _____ **a-D** _____
 (Knee high)

C _____ a-C _____

B' _____ a-B' _____

a-B _____

_____ Y _____

_____ A _____

Foot Length - Left

For open toe _____ cm


For closed toe _____ cm
 (longest toe)

Foot Length - Right


For open toe _____ cm

For closed toe _____ cm
 (longest toe)

Left Foot Slant



Right Foot Slant



Note: medial = inside
 lateral = outside

*Only available with slant foot

Customer Service: **0845 122 3600** Compression Therapy Helpline: **0800 389 8424** Email: **compression.uk@jobst.com** Website: **www.bsnmedical.co.uk** **www.JOBSTUK.co.uk**

By completing this order form, you are confirming that you are aware of your obligation to obtain consent from the patient on the processing of their data for the production of their JOBST® compression garment in accordance with the General Data Protection Regulation. For more information and our Privacy Policy, visit **www.bsnmedical.co.uk**

JOBST® Confidence Ordering Information

Compression Class / Option Code Pieces

JOBST® Confidence Knee High



Basic style

CCL 1 (18-21mmHg)	L1-16-04	1
CCL 2 (23-32mmHg)	L2-17-04	1
CCL 3 (34-46mmHg)	L3-18-04	1

Options

Closed toe	L-A001C	1
SoftFit	L-A016C	1
Slipform	L-A011C	1
Ankle functional zone	L-A010C	1

Compression Class / Option Code Pieces

JOBST® Confidence Thigh High



Basic style

CCL 1 (18-21mmHg)	L1-16-06	1
CCL 2 (23-32mmHg)	L2-17-06	1
CCL 3 (34-46mmHg)	L3-18-06	1

Options

Closed toe	L-A001C	1
5cm silicone band	L-A004C	1
Slipform	L-A011C	1
Ankle functional zone	L-A010C	1
Knee functional zone	L-A015C	1

Options



Ankle Functional Zone

Circumferential functional zone allows greater movement and flexibility at the ankle



Knee Functional Zone

Circumferential functional zone allows greater movement and flexibility at the knee



Slipform

Slipform provides a higher shape to the outside of the garment, at the top, conforming to the limb shape



Deco Line

A decorative seam line can be added to the front of the garment on request



Hallux Valgus (bunion relief zone)*

Creates a looser knitting pattern around the ball of the foot to provide relief to painful bunions



SoftFit Top Band

Soft and secure top band that has been uniquely knitted into the garment top band for a more comfortable secure fit



Initials

Patient initials can be knitted into the garment below the JOBST logo for a personal touch

*Only available with slant foot