



Animal Healthcare Catalogue



About Essity

We are focused on the development of world-class, branded products that offer high quality solutions.

Essity is a global leader in the worldwide healthcare market, specialising in the areas of Wound Care and Fracture Management.

Essity stands for brands such as **Soffban®**, **Tensoplast®**, **Co-Plus®**, **Gypsona®** and **Easifix®** with a long history that you have come to trust.

Whether it is dealing with an acute fracture or simply ensuring that an animal is bandaged effectively, Essity has the products to suit all of your veterinary healthcare needs.

Contact information

For additional support and information just reach out to the animal health team.

Email: animalhealthcare@essity.com



Introducing innovative advanced wound care to veterinary professionals



Cutimed®
Siltec®



Cutimed®
Sorbion®



Leukomed®
Sorbact®



Cutimed®
Gel

Contents

About Essity	2
Contact Information	2
What is Cutimed?	4

Wound Dressings

Exudate Management	
Cutimed Siltec	5
Cutimed Sorbion	6

Infection Management	
Cutimed Sorbact	7
Leukomed Sorbact	8

Wound Bed Management	
Cutimed DebriClean	9
Cutimed Gel	9

Tissue Management	
Cutimed Epiona	10

I.V. Film	
Leukomed I.V. Film	10

Post-Operative Dressings	
Leukomed Control	11
Leukomed	11
Leukomed T	12
Leukomed T plus	12

Casting Materials

Cast Tapes and Splints	
Delta-Cast Elite	15
Delta-Cast Conformable	15
Delta-Cast Soft	16
Gypsona	17
Dynacast Prelude	17

Undercast Padding and Stockinette	
Soffban Natural	18
Soffban Synthetic	18
Delta-Dry	18
Delta Terry-Net C	19

Casting Accessories	
Delta-Net Stockinette	20
Delta Terry-Net	20
Tensoban	20

Wound Contact Layers	
Cutimed Cuticell Classic	12
Cutimed Cuticell Contact	12

Skin Protection	
Cutimed PROTECT Spray	12
Cutimed PROTECT Foam Applicator	12
Cutimed PROTECT Cream	12

Bandages	
Retention	
Easifix	13
Easifix K	13
Elastic Adhesive	
Tensoplast	13
Cohesive	
Co-Plus	14
Easifix Cohesive	14

Wide Area Fixation	
Hypafix	21
Hypafix transparent	21

Tapes	
Leukopor	21
Leukosilk	21
Leukofix	22
Leukoplast Sleek	22
Tensoplast	22
Strappal	22

Swabs	
Leukoplast compress cotton gauze	23
Leukoplast compress non-woven	23
Leukoplast swab abdominal	23

What is Cutimed®?

Developed for the management of chronic and non-healing wounds, the Cutimed wound care range reliably covers all stages of the wound healing process



Cutimed® Debriclean

Mechanical debridement pad for superficial wounds and surrounding skin



Cutimed® Siltec®

Foam dressings for low to high exuding wounds



Cutimed® Gel

Hydrogel to support autolytic debridement for necrotic or unclean wounds



Cutimed® Sorbion®

Superabsorbent dressings for low to excessive exuding wounds



Cutimed® Sorbact®

Safe and effective removal of bacteria to reduce the bioburden in wounds and facilitate the healing process



Cutimed® Epiona

Collagen dressing to help boost chronic wounds

Cutimed® Siltec®

A range of foam dressings designed to manage low to highly exuding wounds

The combination of a silicone wound contact layer and superabsorbent stripes above the foam core mean that Cutimed Siltec is excellent in its exudate management and ensures atraumatic dressing changes.

The range

Secure, comfortable, conformable

Cutimed Siltec offers a wide range of sizes, made of skin-friendly polyurethane foam for reliable exudate management from low to high exudate levels.



Cutimed® Siltec® B
Self-adhesive silicone border, secondary fixation is not necessary




Cutimed® Siltec® Plus
Strong initial tack for that helping hand effect, easy to apply and to re-adjust



Cutimed® Siltec®
Skin protection with a feather tack



Cutimed® Siltec® L
Thin dressing, specifically for lower levels of exudate



	Code	Size	Pack Qty
Cutimed Siltec B (Border) ■ ■ ■ ●●●●-●●●● ▾	73284-00	7.5cm x 7.5cm	1 x 10
	73284-01	12.5cm x 12.5cm	1 x 10
	73284-02	15cm x 15cm	1 x 10
	73284-03	17.5cm x 17.5cm	1 x 5
	73284-04	22.5cm x 22.5cm	1 x 5
	73284-05	10cm x 10cm	1 x 10
	73284-06	10cm x 22.5cm	1 x 10
Cutimed Siltec B Oval ■ ■ ■ ●●●●-●●●● ▾	73284-14	7cm x 10cm (oval)	1 x 10
Cutimed Siltec Plus ■ ■ ■ ●●●●-●●●● ▾	73288-00	5cm x 6cm	1 x 10
	73288-01	10cm x 10cm	1 x 10
	73288-02	10cm x 20cm	1 x 10
	73288-03	15cm x 15cm	1 x 10
	73288-04	20cm x 20cm	1 x 5
Cutimed Siltec ■ ■ ■ ●●●●-●●●● ▾	73285-00	5cm x 6cm	1 x 10
	73285-01	10cm x 10cm	1 x 10
	73285-02	10cm x 20cm	1 x 10
	73285-03	15cm x 15cm	1 x 10
	73285-04	20cm x 20cm	1 x 5
Cutimed Siltec L (Light) ■ ■ ■ ●●●●-●●●● ▾	73283-00	5cm x 6cm	1 x 10
	73283-01	10cm x 10cm	1 x 10
	73283-02	15cm x 15cm	1 x 10

Wound depth

- ▾ Superficial
- ▾ Deep
- ▾ Superficial + deep

Exudate level

- Dry
- Low
- Moderate
- High
- Excessive

Tissue type

- Necrotic
- Sloughy
- Granulating
- Epitheliasing
- Infected

Cutimed® Sorbion®

Superabsorbent dressings

Managing exudate levels of different viscosities can be a challenge when treating wounds. Cutimed offers reliable and innovative solutions for all exudate levels, focusing on keeping a well-balanced moisturisation and avoiding maceration – even under compression.

Cutimed® Sorbion® has demonstrated its efficacy in the treatment of highly exuding wounds. Utilising patented Hydration Response Technology, Cutimed® Sorbion® has set standards in modern exudate management.



Expansion rim

- Broad wound contact thanks to even expansion, even under high saturation

Ultrasonic sealing

- Elimination of adhesives and binders for a particularly well-tolerated wound care - mean less patient allergies

Polymer gelling

- Moist wound environment
- Ensures close contact with the wound bed

Fibre matrix

- Avoids clumping
- Vertical absorption, safe and fast

Fibre matrix with polymers

- High exudate absorption capacity
- Extension of dressing change intervals
- No rewet into the wound
- Reduces the risk of maceration

Polypropylene cover

- Strong tear resistance
- High breathability

Code	Size	Pack Qty
Cutimed Sorbion Sachet Extra		
73234-00	5cm x 5cm	1 x 5
73234-01	7.5cm x 7.5cm	1 x 5
73234-02	10cm x 10cm	1 x 5
73234-03	20cm x 10cm	1 x 5
73234-04	20cm x 20cm	1 x 5
73234-05	30cm x 20cm	1 x 5
Cutimed Sorbion Sachet XL		
73240-01	45cm x 25cm	1 x 5
Cutimed Sorbion Sachet Multi Star		
73237-01	8cm	1 x 10
73237-03	14cm	1 x 10
Cutimed Sorbion Border		
71010-00	10cm x 10cm	1 x 10
71010-02	15cm x 15cm	1 x 10
71010-04	10cm x 23cm	1 x 10
71010-03	20cm x 20cm	1 x 10
71010-07	9 x 11cm (Oval)	1 x 10

Wound depth

- ☐ Superficial
- ☐ Deep
- ☐ Superficial + deep

Exudate level

- Dry
- Low
- Moderate
- High
- Excessive

Tissue type

- ☐ Necrotic
- ☐ Sloughy
- ☐ Granulating
- ☐ Epitheliasing
- ☐ Infected








Cutimed® Sorbact®

Bacteria binding dressings with Sorbact® Technology

Safely remove bacteria -

Killing the bacteria releases endotoxins from dead cells and cell debris is left in the wound. So unlike antimicrobial dressings that actively kill bacteria, Cutimed® Sorbact® removes it, irreversibly binding it to its surface to reduce bioburden and support wound healing. Cutimed® Sorbact® has no known contraindications and a low risk of allergies.

Cutimed® Sorbact® is versatile and effective in a variety of wound stages. It can be used to safely reduce bioburden in the short term, and also in the long-term treatment of already infected wounds.

<p>Cutimed® Sorbact Pad® An absorbent dressing pad covered with a Sorbact Technology coated acetate swab. For wounds with low to moderate exudate levels</p>	
<p>Cutimed® Sorbact Contact® (Formerly Sorbact® Swab) Sorbact Technology coated acetate wound contact layer. Primary wound dressing for either deep or superficial wounds</p>	
<p>Cutimed® Sorbact Ribbon® Made from Sorbact Technology coated cotton. Can be used on superficial and deep wounds such as cavity wounds and fistulas</p>	
<p>Cutimed® Sorbact Hydro® (Formerly Sorbact® Gel) A Cutimed Sorbact swab coated with an amorphous hydrogel for dry, sloughy or low exuding wounds</p>	
<p>Cutimed® Sorbact® Hydroactive A unique combination of a Cutimed Sorbact swab combined with an absorbent hydropolymer gel matrix to maintain a moist wound environment. For low to moderate exudate levels</p>	
<p>Cutimed® Siltec® Sorbact® B Cutimed Siltec Sorbact combines infection management with the excellent fluid handling of Cutimed Siltec. For wounds with low to high exudate levels</p>	
<p>Cutimed® Sorbion® Sorbact® A sterile, bacteria and fungi binding super-absorbent wound dressing. For wounds with high to excessive exudate.</p>	

Code	Size	Pack Qty
Cutimed Sorbact Pad		
72161-16	7cm x 9cm	1 x 5
72162-27	10cm x 10cm	1 x 5
72163-14	10cm x 20cm	1 x 20
Cutimed Sorbact Contact (Formerly Swab)		
72164-33	4cm x 6cm folded (11cm x 16cm unfolded)	1 x 5
72165-62	7cm x 9cm folded (17cm x 27cm unfolded)	1 x 5
Cutimed Sorbact Ribbon		
72166-17	2cm x 50cm	1 x 20
72167-11	5cm x 200cm	1 x 10
Cutimed Sorbact Hydroactive (wound pad)		
72646-03	7cm x 8.5cm (5cm x 6.5cm)	1 x 10
72646-04	14cm x 14cm (10cm x 10cm)	1 x 10
72646-05	19cm x 19cm (15cm x 15cm)	1 x 10
72646-06	24cm x 24cm (20cm x 20cm)	1 x 5
72646-07	14cm x 24cm (10cm x 20cm)	1 x 10
Cutimed Sorbact Hydroactive B (wound pad)		
79933-00	7cm x 8.5cm (5cm x 6.5cm)	1 x 10
79933-02	14cm x 14cm (10cm x 10cm)	1 x 10
79933-03	19cm x 19cm (15cm x 15cm)	1 x 10
79933-04	14cm x 24cm (10cm x 20cm)	1 x 10
Cutimed Siltec Sorbact B (wound pad)		
73251-43	7.5cm x 7.5cm (3.9cm x 3.9cm)	1 x 10
73251-44	12.5cm x 12.5cm (8.2cm x 8.2cm)	1 x 10
73251-45	15cm x 15cm (10cm x 10cm)	1 x 10
73251-46	17.5cm x 17.5cm (12.3cm x 12.3cm)	1 x 5
73251-47	22.5cm x 22.5cm (16.5cm x 16.5cm)	1 x 5
79929-05	Sacrum 17.5cm x 17.5cm	1 x 5
79929-06	Sacrum 23cm x 23cm	1 x 5
Cutimed Sorbact Hydro (Formerly Gel)		
72611-44	7.5cm x 7.5cm	1 x 10
72611-45	7.5cm x 15cm	1 x 10
Cutimed Sorbion Sorbact		
72698-16	10cm x 10cm (8 x 9cm)	1 x 10
72698-17	10cm x 20cm (7.7 x 17.8cm)	1 x 10
72698-18	20cm x 20cm (17.8 x 17.8cm)	1 x 10
72698-19	20cm x 30cm (17.5 x 27.5cm)	1 x 10

Wound depth

- Superficial
- Deep
- Superficial + deep

Exudate level

- Dry
- Low
- Moderate
- High
- Excessive

Tissue type

- Necrotic
- Sloughy
- Granulating
- Epitheliasing
- Infected

Leukomed® Sorbact®

Innovative post-operative dressing to help to reduce wound bacteria colonisation by its physical mode of action. For post-operative and traumatic wounds

Unique bacteria binding Sorbact® Technology

- Effectively binds hydrophobic bacteria and fungi
- Development of bacterial and fungal resistance is not expected

Bacteria proof adhesive film

- Effectively protects against external contamination
- Breathable and waterproof




Code	Size	Pack Qty
76199-28	5cm x 7.2cm	1 x 50
76199-29	8cm x 10cm	1 x 20
76199-30	8cm x 15cm	1 x 20
76199-31	10cm x 20cm	1 x 20
76199-32	10cm x 25cm	1 x 20
76199-33	10cm x 30cm	1 x 20
76199-34	10cm x 35cm	1 x 20



Cutimed® DebricClean

Fast and effective wound debridement

In chronic wounds, devitalised tissue and biofilm forms a physical barrier that inhibits the formation of new tissue. For the wound to heal properly, the devitalised tissue and biofilm must be removed.¹ Cutimed® DebricClean has been designed to support you in the mechanical debridement process and lead to more effective wound healing.



Cutimed® DebricClean

Mechanical debridement pad

- A dual action debridement pad with soft and abrasive fibres
- Effectively removes slough, debris and biofilm
- The hand loop provides easy, comfortable and accurate use

Ideal for: the debridement of superficial wounds and surrounding skin

Code	Size	Pack Qty
76170-00	9cm x 13cm	1 x 5



Cutimed® Gel

Amorphous hydrogel

- A clear, amorphous hydrogel which can be used to debride necrotic and sloughy tissue
- Easy to apply with the sterile nozzle provided, a spatula or directly from the tube
- Adapts to the wound bed

Ideal for: autolytic debridement, sloughy or necrotic tissue, and moisture donation for dry to low exuding wounds


Code	Size	Pack Qty
72610-00	8g	1 x 10
72610-01	15g	1 x 10
72610-02	25g	1 x 10



Cutimed® Epiona

Boosts healing in chronic and stagnating wounds

Natural mode of action - made from 90% native collagen with 10% alginate



Resorbable - no need to remove from the wound

Reduces wound size by up to 40% in 3 weeks¹

- Accelerated wound healing¹
- Safe application^{1,2}
- Can be used with compression therapy¹
- Excellent hemostatic properties^{3,4}
- Broad indication for treatment of all wounds healing by secondary intention, free of necrotic tissue^{1,2}

Code	Size	Pack Qty
73229-00	5cm x 5cm	1 x 10
73229-03	10cm x 10cm	1 x 10



1. Sabo, M, et al., 'A Post-marketing Surveillance Study of Chronic Wounds Treated With a Native Collagen Calcium Alginate Dressing', Ostomy Wound Management, 2018
 2. Cutimed® Epiona - Case Report, Data on File.
 3. Wiegand C., et al., 'A novel native collagen dressing with advantageous properties to promote physiological wound healing', Journal of Wound Care, 2016.
 4. Jesty, J. et al., 'Assessment In Vitro of the Active Hemostatic Properties of Wound Dressings', Journal of Biomedical Materials Research Part B: Applied Biomaterials, 2008

Leukomed® I.V. Film

Sterile transparent dressing for IV sites

- Transparency allows for easy inspection of the puncture site and surrounding area
- Reduces the risk of contamination
- Breathable and can be worn for up to 7 days



Ideal for: securing intravenous catheters and covers and protects cannula sites

Code	Size	Pack Qty
72390-00	6cm x 8cm	1 x 50
72390-03	8.5cm x 11.5cm	1 x 50
72390-04	4.5cm x 4.5cm	1 x 50
72390-05	7cm x 9cm	1 x 50



Leukomed® Control

Transparent post-operative dressing with absorbent hydropolymer pad

- Highly flexible and conforms well to contoured body parts and joints
- Can be worn for up to 14 days depending on the wound condition
- Transparency allows for easy wound monitoring

Ideal for: dry to low exuding acute wounds

Code	Size	Pack Qty
73230-00	5cm x 7cm	1 x 10
73230-01	7cm x 10cm	1 x 10
73230-02	8cm x 15cm	1 x 10
73230-03	10cm x 24cm	1 x 5
73230-05	10cm x 35cm	1 x 5




Leukomed® T

Transparent film dressing

- Allows the skin to breathe
- Reduces the risk of contamination
- Transparency allows for easy wound monitoring

Ideal for: superficial wounds, wounds after surgery and for securing dressings and catheters

Code	Size	Pack Qty
72381-00	7.2cm x 5cm	1 x 50
72381-01	8cm x 10cm	1 x 50
72381-08	10cm x 12.5cm	1 x 50
72381-14	10cm x 25cm	1 x 50
72381-02	11cm x 14cm	1 x 50
72381-09	15cm x 20cm	1 x 50
72381-10	15cm x 25cm	1 x 50



Leukomed®

Non-woven absorbent

- Absorbent and air permeable
- Flexible even on contoured body areas

Ideal for: post-operative care. For short-term use on minor or medium-sized wounds

Code	Size	Pack Qty
72380-00	7.2cm x 5cm	1 x 50
72380-01	8cm x 10cm	1 x 50
72380-02	8cm x 15cm	1 x 50
72380-03	10cm x 20cm	1 x 50
72380-04	10cm x 25cm	1 x 50
72380-05	10cm x 30cm	1 x 50
72380-06	10cm x 35cm	1 x 50



Leukomed® T plus

Transparent film dressing with wound pad

- Ultra thin and conformable film dressing with a low-adherent wound pad
- Sterile barrier protects the wound and supports undisturbed wound healing

Ideal for: dressing surgical incisions, cover and protect superficial cuts, lacerations and minor burns

Code	Size	Pack Qty
72382-00	7.2cm x 5cm	1 x 50
72382-01	8cm x 10cm	1 x 50
72382-02	8cm x 15cm	1 x 50
72382-14	10cm x 20cm	1 x 50
72382-03	10cm x 25cm	1 x 50
72382-04	10cm x 30cm	1 x 50
72382-05	10cm x 35cm	1 x 50



Cutimed® Cuticell® Classic

Paraffin-impregnated gauze dressing

- A low-adherent dressing made from open weave cotton gauze, impregnated with soft paraffin
- Paraffin coating helps to maintain a moist wound environment

Ideal for: exuding wounds

Code	Size	Pack Qty
72538-00	5cm x 5cm - Individually sealed	1 x 5
72538-01	5cm x 5cm - Individually sealed	1 x 50
72538-02	10cm x 10cm - Individually sealed	1 x 10
72538-03	10cm x 10cm - Individually sealed	1 x 100
72538-04	10cm x 40cm - Individually sealed	1 x 10
72538-06	10cm x 7m - Dispenser roll	1 x 1




Cutimed® Cuticell® Contact

Silicone wound contact layer

- Gentle and highly conformable silicone wound contact layer for atraumatic dressing changes

Ideal for: acute and chronic wounds

Code	Size	Pack Qty
72680-00	5cm x 7,5cm	1 x 5
72680-01	7,5cm x 10cm	1 x 5
72680-02	10cm x 18cm	1 x 5
72680-03	15cm x 25cm	1 x 5



Cutimed® PROTECT

Preserving the skin's barrier function

Cutimed® PROTECT Spray

Provides a long-lasting protective barrier against incontinence or exudate, in a spray film format



Cutimed® PROTECT Foam Applicator

Provides a long-lasting protective barrier against incontinence or exudate, in a liquid film format with foam applicator

Cutimed® PROTECT Cream

A moisturising barrier cream that protects skin from the damage and irritation caused by bodily fluids

Code	Size	Pack Qty
Cutimed PROTECT Cream		
72652-00	28g	1 x 1
72652-03	90g	1 x 1
72652-08	2g	1 x 20
Cutimed PROTECT Spray		
72653-00	28ml	1 x 1
Cutimed PROTECT Foam Applicator		
72654-00	1ml	1 x 5
72654-02	3ml	1 x 5

Easifix®

Retention bandage

- Made from a quality viscose and polyamide material to provide a soft, lightweight and porous bandage
- Excellent conformability

Ideal for: fixation of dressings on all body parts especially joints

Code	Size	Pack Qty
71429-10	5cm x 4m stretched	1 x 20
71430-03	7.5cm x 4m stretched	1 x 20
71428-04	10cm x 4m stretched	1 x 20
71429-11	15cm x 4m stretched	1 x 20



Easifix® K

Knitted retention bandage

- Knitted bandage made from viscose and polyamide, giving two-way stretch
- Lightweight, soft and breathable
- Highly elastic and flexible bandage that easily conforms to the most awkward body parts

Ideal for: fixation of dressings on all body parts especially joints

Code	Size	Pack Qty
72617-00	2.5cm x 4m	1 x 20
72617-01	5cm x 4m	1 x 20
72617-02	7.5cm x 4m	1 x 20
72617-03	10cm x 4m	1 x 20
72617-04	15cm x 4m	1 x 20



Tensoplast®

Elastic adhesive bandage

- High quality bandage made from a combination of specially formulated porous adhesive and a superior woven cloth
- Soft, non-fraying edges
- Standard and reverse wound bandage option
- Reverse wound bandage has adhesive part facing outwards


Ideal for: quick, easy and safe application to horses hind legs



Code	Size	Pack Qty
72050-58	2.5cm x 4.5m (plaster)	1 x 12
72050-32	5cm x 4.5m stretched	1 x 12
72050-36	7.5cm x 4.5m stretched	1 x 12
72050-29	10cm x 4.5m stretched	1 x 12
72050-59	5cm x 4.5m stretched (reverse wound)	1 x 60
72050-34	7.5cm x 4.5m stretched (reverse wound)	1 x 48

Tensoplast® Reverse Wound

Elastic adhesive bandage



Tensoplast Reverse Wound bandage is a preferred EAB when bandaging the hind legs of a difficult horse. The adhesive part facing outwards allows for easy roll on application and with little drag effect, making it quicker and therefore safer to apply. This also gives more control with regards to how much tension is applied.



Easifix® Cohesive

Cohesive retention bandage

- Latex-coated cohesive bandage made from viscose and polyester
- Securely sticks to itself for strong fixation of wound dressings and plaster splints
- Only a 30% overlap is required



Ideal for: secure fixation of wound dressings

Code	Size	Pack Qty
71436-14	4cm x 20m stretched	1 x 1
71436-15	6cm x 20m stretched	1 x 1
71436-16	8cm x 20m stretched	1 x 1
71436-17	10cm x 20m stretched	1 x 1
71436-18	12cm x 20m stretched	1 x 1
71436-09	4cm x 4m stretched	1 x 10
71436-10	6cm x 4m stretched	1 x 10
71436-11	8cm x 4m stretched	1 x 10
71436-12	10cm x 4m stretched	1 x 10
71436-13	12cm x 4m stretched	1 x 10

Co-Plus®

Cohesive bandage

- A non-slip, light support bandage suitable for final layer fixation
- Flexible bandage, conforms easily to the most awkward body parts
- Holds in place by sticking to itself

Ideal for: light compression or as secondary fixation




Code	Size	Pack Qty
White		
72101-09	5cm x 6.3m stretched	1 x 36
72101-10	7.5cm x 6.3m stretched	1 x 24
72101-05	10cm x 6.3m stretched	1 x 18
Mixed		
72102-04	5cm x 6.3m stretched	1 x 36
72102-03	7.5cm x 6.3m stretched	1 x 24
72102-02	10cm x 6.3m stretched	1 x 18
Dark blue		
72102-23	10cm x 6.3m stretched	1 x 18
Green		
72102-05	10cm x 6.3m stretched	1 x 18



Delta-Cast® Elite

Non-fibreglass cast tape

- An extremely strong and lightweight cast tape made from polyester substrate
- Extensible yarns which shape easily to body contours
- Soft edges and smooth finish



Ideal for: primary and secondary cast requirements for both small and large animals

Code	Size	Colour	Pack Qty
72279-00	5cm x 3.6m	○ White	1 x 10
72279-01	7.5cm x 3.6m	○ White	1 x 10
72279-02	10cm x 3.6m	○ White	1 x 10
72279-03	12.5cm x 3.6m	○ White	1 x 10

More colours available upon request

To comply to REACH UK guidance, product safety education must be undertaken by the user. To access free training contact: concierge.service@essity.com

Delta-Cast® Conformable

Non-fibreglass cast tape

Delta-Cast® Conformable has set a benchmark in standard and removable, non-fibreglass cast tapes. Delta-Cast® Conformable responds to the daily challenges of fracture management - it's distinctive construction provides rigid and semi-rigid casting options for primary and secondary casting applications where patient fit, comfort and compliance are critical.



Delta-Cast® Conformable is the only UK cast tape with clinical evidence¹ to support its use for focused rigidity casting.

Code	Size	Colour	Pack Qty
72280-00	2.5cm x 1.80m	○ White	1 x 10
72280-01	5cm x 3.6m	○ White	1 x 10
72280-02	7.5cm x 3.6m	○ White	1 x 10
72280-03	10cm x 3.6m	○ White	1 x 10
72280-04	12.5cm x 3.6m	○ White	1 x 10

More colours available upon request

To comply to REACH UK guidance, product safety education must be undertaken by the user. To access free training contact: concierge.service@essity.com

1. Petty, A. C., Wardman, C., (1998). A randomized, controlled comparison of adjustable focused rigidity primary casting technique with standard plaster of Paris/ synthetic casting technique in the management of fractures and other injuries. *Journal of Orthopaedic Nursing*, 2, 95 - 102

Delta-Cast® Soft

Flexible casting tape

Delta-Cast® Soft is designed to provide functional support to fractures and soft tissue injuries. It's semi-rigidity allows for flexibility and mobility, helping to reduce immobilisation related problems such as muscle atrophy and joint stiffness. Improved flexibility and mobility can help promote a faster healing process.

Delta-Cast® Soft is constructed using a combination of high-tech polyester, polypropylene, and elastane yarns, together with a specially developed resin. This technology builds safety, functionality and comfort into a product offering all the advantages only a non-fibreglass casting tape can offer.


Easy removal where the cast can be unpeeled or cut away with cast scissors – no cast saw required.

Ideal for primary and secondary casts requirements, where the cast can be removed without the need of a cast saw.



Code	Size	Colour	Pack Qty
72708-02	2.5cm x 1.8m roll	○ White	1 x 10
72708-12	5cm x 3.6m roll	○ White	1 x 10
72708-22	7.5cm x 3.6m roll	○ White	1 x 10
72708-32	10cm x 3.6m roll	○ White	1 x 10
72708-42	12.5cm x 3.6m roll	○ White	1 x 10
72708-09	5cm x 3.6m roll	🐾 Paw print	1 x 10
72708-19	7.5cm x 3.6m roll	🐾 Paw print	1 x 10
72708-60	5cm x 3.6m roll	❤️ Heart print	1 x 10
72708-61	7.5cm x 3.6m roll	❤️ Heart print	1 x 10

*more colours available upon request



To comply to REACH UK guidance, product safety education must be undertaken by the user. To access free training contact: concierge.service@essity.com



Gypsona®

Premium Plaster of Paris

- Plaster of Paris bandage providing excellent stability and conformability using an interlocked leno weave gauze
- High quality, fast setting plaster of Paris
- Available as a roll or slab



Ideal for: supporting and immobilising fresh fractures and soft tissue injuries

Code	Size	Pack Qty
71985-04	5cm x 2.7m roll	1 x 48
71985-05	7.5cm x 2.7m roll	1 x 48
71985-00	10cm x 2.7m roll	1 x 48
71985-02	15cm x 2.7m roll	1 x 24
71985-03	20cm x 2.7m roll	1 x 24
71988-14	5 ply slab 10cm x 20m	1 x 1
71988-15	5 ply slab 15cm x 20m	1 x 1
71988-16	5 ply slab 20cm x 20m	1 x 1

Dynacast® Prelude Solo

Unpadded precut splinting system

The Dynacast® Prelude Solo unpadded precut is a fibreglass splinting system that comes in a variety of widths and lengths and is stronger and more versatile than traditional plaster of Paris slabs. The splints are easy to apply.

Ideal for treating trauma and post-operative applications.



Code	Size	Pack Qty
72963-00	90cm x 5.5m roll	1 x 10 splints
72963-01	5cm x 2.7m roll	1 x 10 splints
72963-02	7.5cm x 2.7m roll	1 x 10 splints
72963-04	10cm x 2.7m roll	1 x 10 splints
72963-05	15cm x 2.7m roll	1 x 10 splints

To comply to REACH UK guidance, product safety education must be undertaken by the user. To access free training contact: concierge.service@essity.com

Dynacast® Prelude

Synthetic splint system

- A strong, light and durable synthetic splint system incorporating fibreglass, covered with a polypropylene padding
- Quick setting and allows early weight bearing

Ideal for: supporting and immobilising fresh fractures and soft tissue injuries



Code	Size	Pack Qty
71418-00	2.5cm x 4.6m	1 x 1
71420-00	5cm x 4.6m	1 x 1
71421-00	7.5cm x 4.6m	1 x 1
71415-00	10cm x 4.6m	1 x 1
71416-00	12.5cm x 4.6m	1 x 1
71417-00	15cm x 4.6m	1 x 1
71419-00	20cm x 2.3m	1 x 1

To comply to REACH UK guidance, product safety education must be undertaken by the user. To access free training contact: concierge.service@essity.com

Dynacast® Prelude

Strap kits

Designed for use with Dynacast® Prelude splints in the rehabilitation phase, this flexible product enables long-term use without the need for a full bandage, enabling easy splint removal. Each kit contains three latex-free hook and loop straps.

Ideal in fitting splints and removable casts, without the need of a full bandage



Code	Size	Colour	Pack Qty
71415-18	Upper extremity kit	Red	1 x 10
71415-19	Lower extremity kit	Blue	1 x 10

Soffban® Synthetic

Synthetic padding

- Made from 100% polyester
- Autoclavable
- Soft and highly conformable
- Maximum loft retention for excellent cushioning
- Reduced risk of skin maceration and sensitivity
- Tearability for easy application

Ideal for: use under synthetic casts, plaster of Paris and bandages



Code	Size	Pack Qty
71464-01	5cm x 2.7m	1 x 12
71467-01	7.5cm x 2.7m	1 x 12
71486-05	10cm x 2.7m	1 x 12
71465-01	15cm x 2.7m	1 x 12
71466-01	20cm x 2.7m	1 x 6

Soffban® Natural

Natural padding

- Suitable to use where bleeding can occur
- A highly absorbent padding made from 100% viscose rayon
- Excellent conformability
- Easy to tear
- Autoclavable

Ideal for: using as a comfortable undercast layer, use under bandages or under plaster of Paris and synthetic casts



Code	Size	Pack Qty
71472-00	5cm x 2.7m	1 x 12
71473-00	7.5cm x 2.7m	1 x 12
71469-00	10cm x 2.7m	1 x 12
71470-00	15cm x 2.7m	1 x 12
71471-00	20cm x 2.7m	1 x 6

Delta-Dry®

Water-resistant cast padding and stockinette with quick-dry system



- The patented Delta-Dry material offers water-resistance
- The open knit structure of the fibres accelerate drainage after submerging in water
- The body temperature evaporates residual moisture quickly through the airy mesh. This helps the cast and skin dry quickly, skin feels dry in about 90 minutes after daily hygiene
- The polypropylene-nylon-polyester mesh of Delta-Dry products provides a soft, comfortable cushion to the skin
- High breathability keeps skin dry and well aerated helping to reduce skin problems such as maceration, itching and odour

Apply Delta-Dry Using Standard Casting Techniques

- Delta-Dry products are applied the same way as conventional undercast products
- There is no need to learn new application techniques or undergo special training
- The highly conformable and stretchable material makes wrinkle-free application easily and quickly achievable
- Delta-Dry can be removed using standard removal equipment

Ideal to be with synthetic casts to allow animals to be washed down whilst wearing their cast

Code	Size	Pack Qty
Delta-Dry - Water-Resistant Cast Padding		
73443-00	5cm x 2.4m roll	1 x 12
73443-01	7.5cm x 2.4m roll	1 x 12
73443-02	10cm x 2.4m roll	1 x 12
Delta-Dry Stockinette - Water Resistant Stockinette		
74564-00	2.5cm x 10m	1 x 1
74564-01	5.1cm x 10m	1 x 1
74564-02	7.6cm x 10m	1 x 1
74564-03	10.1cm x 10m	1 x 1
74564-04	12.5cm x 10m	1 x 1
Delta-Dry Softliner - Water Resistant Softliner		
72779-01	2.5cm x 10m	1 x 1
72779-02	5.1cm x 10m	1 x 1
72779-03	7.6cm x 10m	1 x 1
72779-04	10.1cm x 10m	1 x 1
72779-05	12.5cm x 10m	1 x 1

Delta Terry-Net® C

Natural terry cloth stockinette

- Cotton, soft stockinette, which provides good cushioning and protection
- Easy to use dispenser
- Autoclavable

Ideal for: underneath padding or synthetic and plaster of Paris casting materials



Code	Size	Pack Qty
72305-00	5cm x 13.7m	1 x 1
72305-01	7.5cm x 13.7m	1 x 1
72305-02	10cm x 13.7m	1 x 1
72305-03	12.5cm x 13.7m	1 x 1
72305-04	15cm x 13.7m	1 x 1



Delta-Net® Stockinette

Synthetic stockinette

Delta-Net® is a soft, comfortable, knitted, synthetic stockinette suitable for use with synthetic or plaster of Paris materials.

Ideal for use as a base layer to protect the skin under the cast



Code	Size	Colour	Pack Qty
72307-05	2.5cm x 22.8m	White	1 x 1
72307-00	5.1cm x 22.8m	White	1 x 1
72307-01	7.6cm x 22.8m	White	1 x 1
72307-02	10.1cm x 22.8m	White	1 x 1
72723-00	2.5cm x 22.8m	Black	1 x 1
72723-01	5.1cm x 22.8m	Black	1 x 1
72723-02	7.6cm x 22.8m	Black	1 x 1
72723-03	10.1cm x 22.8m	Black	1 x 1

Delta Terry-Net®

Adhesive felt liner

Delta Terry-Net® Adhesive Felt Liner is a sheet of 2mm semi compressed felt, used when additional comfort is required.

Ideal in providing additional protection to any bony prominences under a cast



Code	Size	Colour	Pack Qty
72295-05	58cm x 100m roll	White	1 x 1

Tensoban®

Pre-tape underwrap

Tensoban® is a very thin, polyurethane wrap, which can be used under adhesive bandages and tapes to protect the skin in taping and strapping techniques.

Ideal to be used as a pre-tape, applied on top of the padding layer in orthopaedic casting, to provide a smooth, compressed surface before the synthetic casting material is applied



Code	Size	Colour	Pack Qty
71500-08	7cm x 20m roll	White	1 x 12

Hypafix®

Wide area fixation

- Self-adhesive, non-woven wide area dressing fixation
- Easily cut to shape and size required

Ideal for: joints and contoured areas and won't constrict or disturb blood circulation



Code	Size	Pack Qty
71442-03	5cm x 5m	1 x 1
71442-04	10cm x 5m	1 x 1
71443-00	2.5cm x 10m	1 x 24
71443-01	5cm x 10m	1 x 1
71443-02	10cm x 10m	1 x 1
71443-03	15cm x 10m	1 x 1
71443-04	20cm x 10m	1 x 1
71443-05	30cm x 10m	1 x 1

Leukopor®

Breathable medical tape

- A gentle fixation tape made from a non-woven polyester viscose fabric
- Pliable and conforms well to highly contoured areas of the body
- Can be torn easily by hand

Ideal for: sensitive skin to secure dressings, electrodes, probes, tubes, catheters and other measuring devices



Code	Size	Pack Qty
02471-00	1.25cm x 5m	1 x 24
02472-00	2.5cm x 5m	1 x 12
02474-00	5cm x 5m	1 x 6
02453-00	1.25cm x 9.2m	1 x 24
02454-00	2.5cm x 9.2m	1 x 12
02455-00	5cm x 9.2m	1 x 6

Hypafix® transparent

Waterproof wide area fixation

- Thin and conformable film
- Reduces risk of contamination
- Transparency allows for easy inspection of the surrounding wound area

Ideal for: securing small and large wound dressings, devices or instruments



Code	Size	Pack Qty
72378-00	5cm x 10m	1 x 1
72378-01	10cm x 10cm	1 x 1
72378-02	15cm x 10m	1 x 1
72378-03	10cm x 2m	1 x 1

Leukosilk®

Gentle medical tape

- The high strength universal tape made from a fine acetate fabric
- Good initial and permanent adhesion
- Well tolerated by sensitive skin

Ideal for: sensitive skin to secure dressings, tube, catheters, probes and cannulae. For short-term use on intact skin



Code	Size	Pack Qty
72669-00	1.25cm x 4.6m	1 x 24
72669-01	2.5cm x 4.6m	1 x 12
72669-02	5cm x 4.6m	1 x 6

Leukofix®

Transparent medical tape

- Skin-friendly transparent polyethylene film with polyacrylate adhesive
- Provides secure fixation and visibility for wound monitoring
- Water-repellent

Ideal for: securing catheters, probes, cannulae, or fixing occlusive dressings in place



Code	Size	Pack Qty
02121-00	1.25cm x 5m	1 x 24
02122-00	2.5cm x 5m	1 x 12
02124-00	5cm x 5m	1 x 6
02136-00	1.25cm x 9.2m	1 x 24
02137-00	2.5cm x 9.2m	1 x 12
02138-00	5cm x 9.2m	1 x 6

Leukoplast® Sleek

Waterproof, highly adhesive medical tape

- A smooth, thin, plastic film
- Strong, waterproof and reliable
- Latex-free variants available

Ideal for: securing and protecting wound dressings



Code	Size	Pack Qty
72379-01	2.5cm x 3m	1 x 12
72379-02	2.5cm x 5m	1 x 12
72379-03	5cm x 5m	1 x 12
72379-04	7.5cm x 5m	1 x 12
72379-05	10.5cm x 5m	1 x 24

Tensoplast®

Elastic adhesive plaster

- A beige-coloured plain woven fabric
- Easy to unwind
- Can be cut to length

Ideal for: strong and effective, stabilisation of joint structures



Code	Size	Pack Qty
72050-58	2.5cm x 4.5m (plaster)	1 x 12

Strappal®

Rigid strapping tape

- Constructed from a viscose fabric and zinc oxide adhesive
- Hand tearable for quick and secure application

Ideal for: strong and effective stabilisation of joint structures



Code	Size	Pack Qty
71496-15	1.25cm x 10m	1 x 36
71487-00	2.5cm x 10m	1 x 36
71490-00	4cm x 10m	1 x 24
71493-00	5cm x 10m	1 x 18
71494-01	7.5cm x 10m	1 x 6

Leukoplast® compress cotton gauze

Gauze swabs

- Made from high quality gauze
- Available in two colours; green for anaesthetics and white for general use
- Available sterile and non-sterile

Ideal for: wound cleansing, skin preparation, general cleansing or as a primary or secondary dressing

Code	Size	Sterile	X-ray	Type	Ply	Colour	Pack Qty
71237-00	5cm x 5cm	✓	✗	T17	8	White	50 (25 x 2)
71237-08	5cm x 5cm	✓	✗	T17	12	White	50 (25 x 2)
71238-00	7.5cm x 7.5cm	✓	✗	T13 Heavy	8	White	240 (48 x 5)
71237-02	7.5cm x 7.5cm	✓	✗	T17	8	White	50 (25 x 2)
71237-11	7.5cm x 7.5cm	✓	✗	T17	12	White	50 (25 x 2)
71237-04	10cm x 10cm	✓	✗	T17	8	White	50 (25 x 2)
71237-14	10cm x 10cm	✓	✗	T17	12	White	50 (25 x 2)
71233-00	10cm x 10cm	✓	✓	T17	12	White	200 (20 x 10)
71233-07	10cm x 10cm	✓	✓	T17	12	White	140 (14 x 10)
71235-00	5cm x 5cm	✗	✗	T13	8	White	100
71235-05	5cm x 5cm	✗	✗	T13	12	White	100
71235-02	7.5cm x 7.5cm	✗	✗	T13	8	White	100
71235-01	7.5cm x 7.5cm	✗	✗	T13	8	Light green	100
71235-08	7.5cm x 7.5cm	✗	✗	T13	12	White	100
71235-04	10cm x 10cm	✗	✗	T13	8	White	100
71235-06	10cm x 10cm	✗	✗	T13	12	White	100
71235-07	10cm x 10cm	✗	✗	T13	12	Light green	100



Leukoplast® compress non-woven

Non-woven swabs

- Breathable, soft and absorbent swabs
- Available sterile or non-sterile

Ideal for: cleaning wounds, skin preparation or as a secondary dressing

Code	Size	Sterile	X-ray	Type	Ply	Colour	Pack Qty
71288-03	5cm x 5cm	✓	✗	30gr	4	White	250 (50 x 5)
71288-00	5cm x 5cm	✓	✗	30gr	4	White	100 (50 x 2)
71288-01	5cm x 5cm	✓	✗	30gr	6	White	75 (75 x 1)
71288-06	7.5cm x 7.5cm	✓	✗	30gr	4	White	250 (50 x 5)
71288-04	7.5cm x 7.5cm	✓	✗	30gr	4	White	100 (50 x 2)
71288-05	7.5cm x 7.5cm	✓	✗	30gr	4	White	12 (12 x 1)
71288-13	10cm x 10cm	✓	✗	30gr	4	White	250 (50 x 5)
71288-09	10cm x 10cm	✓	✗	30gr	4	White	100 (50 x 2)
71288-12	10cm x 10cm	✓	✗	30gr	6	White	75 (75 x 1)
71288-11	10cm x 10cm	✓	✗	30gr	4	White	6 (6 x 1)
71287-00	5cm x 5cm	✗	✗	30gr	4	White	100
71287-01	7.5cm x 7.5cm	✗	✗	30gr	4	White	100
71287-02	7.5cm x 7.5cm	✗	✗	30gr	6	White	50
71287-03	10cm x 10cm	✗	✗	30gr	4	White	100
71287-04	10cm x 10cm	✗	✗	30gr	6	White	50



Leukoplast® swab abdominal

Sterile X-Ray detectable swabs

- Made from high quality cotton gauze
- Available in two colours, green and white

Ideal for: surgical procedures

Code	Size	Sterile	X-ray	Ply	Colour	Pack Qty
71231-09	20cm x 30cm	✓	✓	2	Green	70 (14 x 5)
71232-11	20cm x 30cm	✓	✓	2	White	70 (14 x 5)
71231-05	45cm x 45cm	✓	✓	2	Green	30 (6 x 5)
71232-05	45cm x 45cm	✓	✓	2	White	30 (6 x 5)



Education by Essity

The Essity Animal Health team are excited to be able to conduct wound care education and product training with you via online meetings.

Training topics covered include;

- General wound care products
- Debridement
- Exudate management
- Post-operative wound protection

We can tailor the training to suit your needs.

Book a training session with us and you will receive a free advanced wound care sample pack which will contain items specific to your training.

Contact the team: animalhealthcare@essity.com

Please contact the team for more information

Tel: **01482 670100**

Email: animalhealthcare@essity.com